
m/r IGF-I - ELISA

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Enzyme Immunoassay for the Quantitative Determination of
Mouse-/Rat- Insulin-like Growth Factor-I
(IGFBP-blocked)
English

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For Research Use Only.
Not for use in diagnostic procedures.

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according to DIN EN 980 and EDMA recommendations Standard News 6 2001

Expiry date/ Verfallsdatum/ Date de péremption/ Data di scadenza/ Fecha de caducidad/ Data de validade/ Uiterste gebruiksdatum/ Udløbsdato/ Bäst före-datum/ Termin ważności/ Lejárati idő/ Čas expirácie/ Doba expirace/ Срок на годност/ Αεγυμίσκυρπäv/ Ημερομηνία λήξης/ Data de expirare/ Rok uporabe/ Viimeinen käyttöpäivä

Consider instructions for use/ Bitte Gebrauchsanweisung beachten/ Consultez la notice d'utilisation/ Consultare le istruzioni per l'uso/ Consulte las instrucciones de uso/ Respeitar as instruções de utilização, / A.u.b de gebruiksaanwijzing volgen/ Se brugsanvisningen/ Läs anvisningarna före användning/ Proszę przeczytać instrukcję obsługi/ Vegye figyelembe a használati utasításban foglaltakat/ Postupujte podľa pokynov na použitie/ Dodržujte návod k použití/ Моля, спазвайте инструкцията за употреба/ Palun järgige kasutusjuhendit./ Λάβετε υπόψη σας τις οδηγίες χρήσης/ Vá rugám sã respectați instrucțiunile de utilizare/ Upošteвайте navodila za uporabo/ Lue käyttöohje huolellisesti!

Lot-Batch Number/ Charge-Chargennummer/ Lot-Code du lot/ Lotto-Numero di lotto/ Lote-Código de lote/ Lote-Código do Lote/Lot-Partijnummer/ Lot-Batchcode/ Partisatskod/ Numer serii/ Tétel-sarzs szám/ Číslo šarže/ Číslo šarže/Партиден номер/Partii – partii number/ Παρτίδα-αριθμός παρτίδας/ Lot-număr lot/ Številka serije/ erä

Manufactured by/ Hergestellt von/Fabriqué par/ Prodotto da/Fabricado por/ Fabricado por/ Vervaardigd doo/Fabrikation af /Tillverkad av/ Wyprodukowane przez/ Gyártotta / Vyrobené/ Vyrobeno v/ Производител/ Тootja/Κατασκευάζεται από/ Produs de/Proizvajalec/ Valmistaja

Catalogue Number/ Bestellnummer/ Numéro de référence/Numero di riferimento/Número de referencia/ Número de Referência/ Referentienummer/Referencennummer /Beställningsnummer/ Numer katalogowy/ Rendelési szám/Katalógovné číslo/ Objednací číslo/Каталожен номер/Tellimisnumber/ Αρ. παραγγελίας/Număr comandă/ Številka naročila/ viite tai tilausnumero

Store at between/ Lagerung bei zwischen/ Conserver à entre/ Conservare a tra/ Conservar a temp. Entre/ Armazenar entre/ Bewaar bij tussen/ Opbevaars mellem/ Förvaras vid/ Przechowywać w/ Tárolási tartomány/ Skladujte v rozsahu / Skladujte v rozmezí/ Температурно ограничение/ Säiläta temperatuuridel/ Φύλαξη σε θερμοκρασία/ Depozitare între/ Skladištenje med/ Säilytys x-y Celsiusasteen lämpötilassa

Contains sufficient for x tests/ Inhalt ausreichend für x Tests/ Contient suffisant pour x tests/ Contenido suficiente para x tests/ Bevat voldoende voor x bepalingen/ Indeholder tilstrækkeligt til x prøver/Innehållet räcker till x analyser/ Zawartość na x testów/ Tartalma x teszt elvégzésére elegendő/ Obsahuje materiál pre x testov / Obsah dostacuje pro x testů/ Съдържание достатъчно за x тестове/ Sisust jätkub x katse jaoks/ Το περιεχόμενο επαρκεί για x δοκιμές/ Conținut suficient pentru x teste/Vsebina zadostuje za x preizkusov/ Sisältö riittää x testille

Keep away from sunlight/ Nicht dem Sonnenlicht aussetzen/ Conserver à l'abri de la lumière/ Conservare al riparo della luce solare/ No exponer a la luz solar/ Proteger da luz solar/ Niet aan zonlicht blootstellen/ Må ikke udsættes for sollys/ Utsätt inte för solljus/ Nie wystawiać na słońce/ Napfénytől távol tartandó/ Nevystavovat slnečnému svetlu/ Nevystavovat slunečnímu světlu/ Да се предпазва от слънчева светлина/ Kaitsta otsese päikesekiirguse eest/ Κρατήστε το μακριά από την ηλιακή ακτινοβολία/ Țineți departe de lumina soarelui/ Ne izpostavljajte sončni svetlobi/ suojaa auringonvalolta

Incubation time/ Inkubationszeit/ Temps d'incubation/ Tempo d'incubazione/ Tiempo de incubación/ Tempo de incubação/ incubatietijd/Inkubationstid/ inkubationstid/ Czas inkubacji/ Inkubációs idő/ Inkubačná lehota/ Inkubační doba/ Инкубационен период/ Inkubatsiooniaeg/ Χρόνος επώασης/ Timp de incubare/Inkubacijska doba/ inkubaatioaika

incubate at / Inkubation bei/ Incuber à/Incubare a/incubar a/Incubar a/incubatietemperatuur/Inkubation ved/inkubation vid/Inkubacja przy/Inkubáció hőmérséklete/Inkubácia pri/Inkubace při/Инкубира се при/Inkubatsioon temperatuuril/Επώαση στους/Incubare la/Inkubacija pri/ inkubaatiolämpötila

Shaking/ Schütteln/ Mélanger/ Agitare/ Agitar/ Agitação/ Schudden/ Ryster/ Skaka/ Wstrząsanie/ Rázás/ Pretrepat/ Protřepat/ Разклащане/ Raputada/ Ανακινήστε/ Vibrare/ Stresite/ Sekoita

Microtiterplate/ Mikrotiterplatte/ plaque de microtitrage/ Piastra di microtitolazione/ Placa de microtitulación/ Placa de Microtitulação/ microtiterplaat/ Mikrotiterplade/ mikrotiterplatta/ microtiterplaat/ Płytko microtiter/ Mikrotiter lap/ Mikrotitračná podložka/ Mikrotitrační podložka/ Микротитърна плака/ Mikrotiiterplaat/ Τρυβλίο μικροπιλοδότησης/ Microplacă/ Mikrotitrská plošča/ Mikrotitrauslevy

Reconstitute in/ Rekonstituieren in/ Reconstituer dans/ Ricostituire nel/ Reconstituer en/ Reconstituer em/ reconstituieren in/ Rekonstituér i/ rekonstituera/ Rekonstituować w/ Helyreállítás/ Znovu připravit za/ Znovu připravit za/ Разтваряне в/ Moodustada uuesti / Ανασυστήστε σε/ Reconstituire în/ Predelava v /rekonstituoi

Sample/ Probe /Echantillon/ campione/ Muestra/ Amostra/ monster/ Prøve/ prov/ Próbka/ Minta/ Vzorka/ Vzorek/ Проба/ Proov/ Δείγμα/ Probă/Vzorec/Näyte

Antibody Conjugate ; Antikörperkonjugat, anticorps conjugué, Coniugato di anticorpo, Conjugado de anticuerpos, Conjugado Anticorpo, antilichaamconjugaat, Antistoffer-konjugat, antikroppskonjugat, Konjugat antyciał, Antitest páros, Prottilátkový konjugát, Prottilátkový konjugát, Антигено конюгат, Antikehad konjugaat Σύμπλοκο αντισώματος, Compuși din anticorpi, Antitelesa konjugat, Vasta-aine konjugaat

Enzyme Conjugate, Enzym Konjugat, conjugué enzymatique, Coniugato di enzima, Conjugado de enzimas, Conjugado Enzima, enzymconjugaat, enzym-konjugat, enzymkonjugat, Konjugat enzymów, Enzim páros, enzymatický konjugát, enzymatický konjugát, ензим конюгат, Ensüümi konjugaat, Σύμπλοκο –ενζύμου, Compuși din enzime, Encima Konjugat, enzimi-konjugaat

Dilute in Buffer X/ Verdünnen in Puffer X/ Diluer dans le tampon X/ Diluire nel tampone X/ Diluir en tampón X/ Diluir no Tampão X/ verdunnen in buffer X/ Fortyndes i buffer X/ späd i buffert X/ Rozcieńczenie w buforze X/ Hígítás X pufferben/ Riedit' v pufrí X/ Ředit v pufru X/ Разреждане в буфер X/ Lahjendada puhvis X/ Αραιώστε σε ρυθμιστικό διάλυμα X/ Diluați în tamponul X/ Razredčiti v pufru X/ laimennetaan x puskuriin

Standard X /Standard X/ Etalon X/ Standard X/ Estándar X/ Standard X/ standaard X/ Standard X/ standard X/ Standard X/ Standard X/ Štandard X/ Standard X/ Стандарт X/ Standard X/ Πρότυπο X/ Standard X/

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For Research Use Only!

For professional use only!

CAUTION: Not for human or animal therapeutic or diagnostic use.

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Package Insert English

Table of Content

Intended Use	16
Reagents Provided	17
Materials Required but not provided	18
WARNINGS AND PRECAUTIONS	18
PRINCIPLE	19
Specimen	19
Storage of the samples	20
Sample Preparation	20
Technical Notes	20
Standards and Controls	21
Washing Buffer	21
Microtiterplate	21
Assay Procedure	22
Calculation of Results	23
Establishing the Standard Curve	23
Performance Characteristics	24
Standards	24
Sensitivity	24
Interassay Variability	24
Intraassay Variability	24
Linearity	25
LITREATURE / LITERATUR	26
KURZANLEITUNG – Mediagnost IGF-I-ELISA E25	28
SUMMARY – MEDIAGNOST IGF-I ELISA E25	30
REF E25 International Test Description	32

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
PACKAGE INSERT ENGLISH

Mediagnost Mouse/Rat IGF-I ELISA E25

- is suited for IGF-I determination in serum and plasma of **mice and rats**
- high sensitive: **0.029 ng/ml** analytical sensitivity; required sample volume is very small
- is **fast**: incubation time a total of 2 hours
- Single Standards with **0.5, 2.5, 6, 12, 18 ng/ml** recombinant IGF-I are provided in the Kit
- 2 Control Sera are provided for quality control
- uses **high affinity antibodies** against m/r IGF-I
- Microtiter plates are separately breakapart

Intended Use

Measurement of IGF-I in mouse / rat serum and plasma.

INTRODUCTION

Beside different cell culture models and studies with human patients, mice and rats are suitable model organisms for basic research and pre-clinical studies. Thus, we developed this test system as a tool for IGF-I measurements in mice and rat for usage in research and pre-clinical studies. Even if the comparability of mice and humans is limited we offer some background information on the *human* IGF-I system in the following section:

Insulin-like growth factors (IGF) I and II play a pivotal role in regulating the proliferation and differentiation of many cell types (1-3). IGF-I is identical with Somatomedin C (Sm-C) (4) and has a molecular weight of 7649 daltons (5). Its major regulators are growth hormone (GH) and nutrition (6). In contrast to many other peptide hormones, IGFs are avidly bound to specific binding proteins (IGFBP). The seven IGFBPs which are known at present (7,8,22) either bind IGF-I and IGF-II with similar affinities or show a preference for IGF-II (9,10).

A major problem of IGF-I measurement results from the interference of IGFBPs in the assay. Direct determinations in untreated serum samples (11) give false values because of the extremely slow dissociation of the IGF-I/IGFBP-3 complexes during the assay incubation. Depending on the ratio IGF-I to IGFBP in the sample interference comes up (see example Figure 1):


Figure 1. Interference of IGFBP in IGF-I measurements. Known concentrations of IGF-I were assayed in the presence of 0.5 ng (left) or 5 ng (right) hIGFBP-3 by a conventional (□) and by the IGFBP-blocked assay (*).

Therefore, various techniques were applied to physically separate IGF-I from its binding proteins before measurement, including (a) size exclusion chromatography under acidic conditions, (b) solid-phase extraction and (c) acid-ethanol extraction (2,12,13). These techniques, however, are either inconvenient or time-consuming or give incomplete and not-reproducible recoveries. The most widely used method is the acid-ethanol extraction (13,14) with a recovery of only 70-80 % of IGFBP-bound IGF-I as a result of co-precipitation. The absolute results of such an extraction are therefore false low (15). The extraction removes the IGFBPs only insufficiently and leads to reduction in sensitivity of the assay due to pre-dilution of the samples by the extraction procedure. Furthermore, the remaining IGFBP may still interfere in the assay. In addition, the acid-ethanol extraction is ineffective in specimens other than serum or plasma (e.g. cell culture media), in which determination of IGF-I is already difficult enough due to the fact that IGFBPs are frequently present at large excess.

To avoid these difficulties, an uncomplicated assay was developed, in which special sample preparation is not required before measurement.

Reagents Provided

1)	MTP	Microtiter plate , ready for use: Microtiter plate with 96 wells, divided up in 12 stripes à 8 wells (separately breakapart) coated with anti-mouse/rat IGF-I antibody.
2)	CAL	Standards A-E , lyophilised, contain recombinant IGF-I. Standard values are between 0.5 - 18 ng/ml (0.5, 2.5, 6, 12 und 18 ng/ml) IGF-I and have to be reconstituted in 1 ml (each) Sample Buffer PP . 50 µl per well are used in the assay. If the standards are required for more than one assay process we recommend to store the reconstituted Standards frozen at -20°C. Attention: Standards should be thawed only once – where required please store aliquoted in adequate volumes.
3)	DILU	Sample buffer PP, 125 ml , ready for use, please use for the reconstitution of Standards A-E, Control Sera KS1 & KS2 and for the serum dilution.
4)	Control	Control Sera KS1 & KS2 , each 500 µl lyophilised: Sera have to be reconstituted in 500 µl Sample Buffer PP . The IGF-I target values and the respective ranges are given on the vial labels. The dilution of the Control Sera KS 1&2 should be according to the dilution of the respective samples, the target values should be obtained by multiplication with the respective dilution factor .
6)	Ab	Antibody Conjugate AK, 7 ml , ready for use, contains the biotinylated anti-IGF-I antibody. Use 50 µl for each well in the assay.
7)	CONJ	Enzyme Conjugate EK, 12 ml , ready for use, contains horseradish-peroxidase conjugate to streptavidin, use 100 µl for each well in the assay.
6)	WASHBUF 20x	Washing Buffer (WP), 50 ml, 20 X concentrated solution. Washing Buffer (WP) has to be diluted 1:20 with distilled or demineralised water before use (e.g. add the complete contents of the flask (50 ml) into a graduated flask and fill up with A.dest. to 1000 ml). Attention: After dilution the Washing Buffer is only 4 weeks stable, dilute only according to requirements.
7)	SUBST	Substrate (S), 12 ml , ready for use, horseradish-peroxidase-(HRP)-substrate, stabilised H ₂ O ₂ Tetramethylbencidine.
8)	H ₂ SO ₄	Stopping Solution (SL), 12 ml , ready for use, 0.2 M sulphuric acid, Caution acid!
9)		Sealing tape for covering of the microtiter plate, 2 x, adhesive.

Materials Required but not provided

Precision pipettes and multichannel pipettes with disposable plastic tips
 Vortex-mixer
 Microtiter plate shaker (350 rpm)
 Microtiter plate washer (recommended)
 Micro plate reader ("ELISA-Reader") with filter for 450 and 620nm (or ≥ 590 nm)
 Polyethylene PE/Polypropylene PP tubes for dilution of samples

WARNINGS AND PRECAUTIONS**For in vitro use only. For professional use only.**

Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.

Before use, all kit components should be brought **to room temperature at 20 - 25°C**. Precipitates in buffers should be dissolved before use by thorough mixing and warming.

Do not mix reagents of different lots. Do not use expired reagents.

The microplate contains snap-off strips. Unused wells must be stored at 2 - 8°C in the sealed foil pouch and used in the frame provided.

Caution: This kit contains material of human and/or animal origin.

No known test methods can offer total assurance of the absence of infectious agents; therefore all components and specimens should be treated as potentially infectious.

Following components contain < 0.01% 2-Methyl-4-isothiazolin-3-one solution as preservative: **A-E, AK, EK, PP**

R34 Irritating to eyes and skin
 R43 Sensibilisation through skin contact possible
 S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
 S36/37 Wear suitable protective clothing and gloves
 S45 In case of accident or if you feel unwell seek medical advice

Following components contain < 0.01% (w/w) 5-chloro-2-methyl 2H isothiazol-3-one and 2-methyl-2H-isothiazol-3-one as preservative: **A-E, AK, EK, PP, WP**

R36/38 Irritating to eyes and skin
 R43 Sensibilisation through skin contact possible
 S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
 S28.1 After contact with skin, wash immediately with plenty of water

Stop solution contains 0.2 M Sulfuric Acid (H₂SO₄)

R36/38 Irritating to eyes and skin
 S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
 S28.1 After contact with skin, wash immediately with plenty of water
 S36/37 Wear suitable protective clothing and gloves.

Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step. Use separate pipette tips for each sample, control and reagent to avoid cross contamination. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.

TMB-Substrate (S) contains 3,3',5,5' Tetramethylbenzidine. Store and Incubate in the dark.

R20/21/R22 Harmful by inhalation, in contact with skin and if swallowed
 R36/37/38 Irritating to eyes, respiratory system and skin
 S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
 S28.1 After contact with skin, wash immediately with plenty of water
 S36/37 Wear suitable protective clothing and gloves

General first aid procedures:

Skin contact: Wash affected area thoroughly with water. Discard contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: If swallowed, wash out mouth thoroughly with water. Immediately see a physician.

Do not eat, drink or smoke in these areas.

Never pipette the materials with the mouth.

Spilled material must be wiped off immediately and should become disinfected. Clean contaminated areas and equipment with a suitable detergent.

PRINCIPLE

The Mediagnost ELISA for m/r IGF-I E25 is a so-called Sandwich-Assay. It utilizes two specific and high affinity antibodies for this protein. The IGF-I in the sample binds to the immobilized first antibody on the microtiter plate, the biotinylated and Streptavidin-Peroxidase conjugated second specific anti-IGF-I-Antibody binds in turn to the immobilized IGF-I. In the closing substrate reaction the turn of the colour will be high specific catalysed, quantitatively depending on the IGF-I-level of the samples.


Figure 2: Principle of the IGFBP blocked IGF-I ELISA

In order to dissociate IGF-I from the IGFBPs, the samples must be diluted in an acidic buffer (Figure 2). The diluted samples are then pipetted into the wells, by this the pH-value will be neutralized. After neutralization of the samples, the excess IGF-II occupies the IGF-binding sites of the binding proteins, thus allowing the measurement of resulting free IGF-I. With this method, the IGFBPs are not removed, but their function and therefore their interference in the assay is neutralized. Due to the extremely low cross-reactivity of the IGF-I antibody with IGF-II, the excess of IGF-II does not disturb the interaction with IGF-I.

The test runs like a conventional ELISA using a Streptavidin-Peroxidase-Enzyme Conjugate.

Specimen

Serum samples as well as Heparin-, EDTA- and Citrat-Plasma samples are suited. Possible dilution of the sample by the anticoagulant must be considered.

Influence of Heparin (30 IE/mL), EDTA (6,8 mM) and NaCitrat (0,015 M) on the measurement of IGF-I has been investigated in recovery experiments. Buffer solution was enriched with recombinant IGF-I and the above mentioned substances. No significant influence on the recovery of IGF-I was detected, on average the recovery of recombinant material in comparison to enriched PBS was 108%.

Cell culture medium is suitable as sample matrix after predilution of 1:2 with Sample Buffer PP.

Haemolytic reactions have to be avoided. The blood has to be allowed to clot and after complete clotting, serum is separated by centrifugation.

Storage of the samples

Storage at RT max. 2 days

Storage at -20°C max. 2 years

More than 2 freeze/thaw cycles are not recommended. IGF-I in samples was found to be unstable under repeated freeze/thaw cycles, measured values declined respectively.

Sample Preparation

Samples have to be diluted at least 1:10 in Sample Buffer (PP).

A serum dilution of 1:100 is in general suitable. However, the IGF-I levels can vary individually significant, we would therefore recommend to check this and adjust the dilution respectively.

Technical Notes

The assay has to be conducted strictly according the test protocol herein.

Reagents with different lot numbers cannot be mixed. The microtiterplate and reagents are stable until the indicated expiry if stored unopened and protected from sunlight at 2 – 8°C.

Bring all reagents to room temperature (20 - 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming.

Incubation at room temperature means: 20-25°C

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtitre plate shaker. We are recommending 350 rpm. Due to certain technical differences deviations may occur, in case the rotation frequency must become adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/or false values, excessive shaking may result in high optical densities and/or false values.

Proper washing is of basic importance for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided washing buffer diluted to usage concentration. Washing volume per washing cycle and well must be 300 µl at least.

The danger of handling with potentially infectious material must be taken into account.

When using an automatic microtitre plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non fuzzy absorbent tissue.

Manual washing is an adequate alternative option. Washing Buffer may be dispensed via a multistep device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamically swinging out the microtitre plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non fuzzy absorbent tissue.

Standards and Controls

For the reconstitution of the lyophilised components (Standards A - E and Control Sera KS1 & KS2) the kit Sample Buffer PP has to be used. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam!) with a Vortex mixer.

The reconstituted standard and controls can be stored for 2 months at -20°C . Repeated freeze/thaw cycles have to be avoided.

Washing Buffer

The required volume of washing buffer is prepared by 1:20 dilution of the provided 20fold concentrate with deionised water. The diluted Washing Buffer is stable for 4 weeks at $2-8^{\circ}\text{C}$. It has to be at room temperature for usage!

Microtiter plate

Unused microtiter plate stripes have to be stored airtight together with the desiccant bag at $2-8^{\circ}\text{C}$. The labelled expiry is not influenced in case of proper storage.

Assay Procedure

All determinations (Standards, Control Sera KS1 & KS2 and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

When performing the assay, the Standards, Control Sera and the samples should be pipette as fast as possible (e.g., <15 minutes).

All incubations have to be conducted at room temperature (20-25°C)

To avoid distortions due to differences in incubation times, Antibody (AK) and Enzyme Conjugate (EK) as well as the following Substrate Solution S should be added to the plate in the same order and in the same time interval as the samples. Stop Solution SL should be added to the plate in the same order as the Substrate Solution.

- 1) Add **50 µl Antibody Conjugate AK** in **all** wells used
- 2) Add **50 µl Sample Buffer PP** into the first two wells (These wells serve as blanks). Subsequently add 50 µl Standard or 50 µl of diluted Control Sera or diluted samples into wells designated for standards, controls or samples.
- 3) Cover the wells with sealing tape and incubate the plate for **1 hour** at **room temperature** (shake at 350 rpm)
- 4) After incubation, aspirate the contents of the wells and wash the wells **5 times 300 µl Washing Buffer WP / well**. The washing buffer should incubate for at least for 15 seconds/cycle. Tap the plate firmly against several layers of folded paper towels to remove excess washing buffer.
- 5) Following the last washing step pipette **100 µl** of the **Enzyme Conjugate EK** in each well.
- 6) Cover the wells with sealing tape and incubate the plate for **0.5 hour** at **room temperature** (shake 350 rpm).
- 7) After incubation wash the wells **5 times** with Washing Buffer as described in step 4
- 8) Pipette **100 µl of the TMB Substrate** Solution in each well.
- 9) Incubate the plate for **30 minutes in the dark at room temperature (20 - 25°C)**.
- 10) Stop the reaction by adding **100 µl of Stopping Solution**.
- 11) Measure the colour reaction within 30 minutes at 450nm (reference filter ≥590 nm).

Calculation of Results

Establishing the Standard Curve

For the evaluation of the assay it is preconditioned that the absorbance values of the blank should be below 0.25, these of standard E should be above 1.0.

Samples, which yield higher absorbance values than Standard E are beyond the standard curve, for reliable determinations these samples should be tested anew with a higher dilution.

Standards are provided in the following concentrations (use the concentration unit as preferred):

Standard	A	B	C	D	E
ng/ml	0.5	2.5	6	12	18

- 1) Calculate the **mean absorbance** value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the blank from the mean absorbances of all other values.
- 3) Plot the standard concentrations on the x-axis versus the mean value of the absorbance of the standards on the y-axis.
- 4) Recommendation: Calculation of the standard curve should be done by using a computer program. **A higher-grade polynomial, or four parametric logistic (4-PL) curve fit or non-linear regression** are usually suitable for the evaluation (as might be spline or point-to-point alignment in individual cases).
- 5) The IGF-I concentration of the diluted sample or the diluted control sera KS1&2 in ng/ml (or µg/ml according the chosen unit for the standards) is calculated in this way, the IGF-I concentration of the **undiluted sample** and of KS1 & KS2 is calculated **by multiplication** with the respective dilution factor.

The exemplary shown standard curve in Fig.3 **cannot be used** for calculation of your test results. You have to establish a standard curve for each test you conduct!

Exemplary calculation of the IGF-I concentration of an diluted sample:

Measured extinction of your sample	0.70
Measured extinction of the blank	0.02

Your **measurement programm** will calculate the IGF-I concentration of the diluted sample automatically by using the difference of sample and blank for the calculation. You only have to determine the most suitable curve fit (here: polynomial 3rd degree).

In this exemplary case the following equation is solved by the programm to calculate the IGF-I concentration in the sample:

$$y = -0.031924 + 0.2606x - 0.010994x^2 + 0.0019791x^3$$

$$3.0921 = x$$

if the dilution factor (**1:100**) is taken into account the IGF-I concentration of the undiluted sample is 309.21 ng/mL


Fig. 3: Exemplary Standard Curve with a polynomial 3rd degree as curve fit.

Performance Characteristics

Standards

The Standards of the ELISA E25 are prepared from **recombinant IGF-I** in concentrations of **0.5, 2.5, 6, 12, 18 ng/mL**.

Sensitivity

The analytical sensitivity of the ELISA E25 yields **< 0.029 ng/mL**

Interassay Variability

	Sample 1	Sample 2	Sample 3
Mean	471	744	347
SD	24.81	41.8	23.3
VC%	5.3	5.6	6.72
n	10	10	16

Intraassay Variability

	Sample 1	Sample 2	Sample 3
Mean [µg/mL]	165	794	514
SD	8,0	47	22
VC%	4,9	5,9	4,2
n	10	10	10

Linearity

ng/mL	Probe 1	Probe 2
1:50	261	888
1:75	332	909
1:100	315	980
1:200	281	960
1:400	290	994
1:600	298	997
1:800	325	983


Method Comparison

Figure 4: Method comparison of the Mediatechn E25 m/r IGF-I ELISA and the IDS m/r IGF-I ELISA.

LITREATURE / LITERATUR

- 1) Baxter RC. 1986 The somatomedins: insulin-like growth factors. *Adv Clin Chem.* 25:49-115
- 2) Daughaday WH, Rotwein P. 1989 Insulin-like growth factors I and II. Peptide, messenger ribonucleic acid and gene structures, serum, and tissue concentrations. *Endocr Rev.* 10:68-91
- 3) Spencer EM (Ed.) 1991 *Modern Concepts of Insulin-Like Growth Factors*. New York: Elsevier.
- 4) Klapper DG, Svoboda ME, Van Wyk JJ. 1983 Sequence analysis of somatomedin-C: confirmation of identity with insulin-like growth factor-I. *Endocrinology.* 112:2215-2217.
- 5) Rinderknecht E, Humbel RE. 1978 The amino acid sequence of human insulin-like growth factor I and its structural homology with proinsulin. *J Biol Chem.* 253:2769-2276.
- 6) Clemmons DR, Van Wyk JJ. 1984 Factors controlling blood concentration of somatomedin C. *Clin Endocrinol Metab.* 13:113-143.
- 7) Ballard J, Baxter R, Binoux M, et al. 1989 On the nomenclature of the IGF binding proteins. *Acta Endocrinol (Copenh).* 121:751-752.
- 8) Drop SLS. 1992 Report on the nomenclature of the IGF binding proteins. *J. Clin Endocrinol Metab.* 74:1215-1216.
- 9) Martin JL, Baxter RC. 1986 Insulin-like growth factor binding protein from human plasma. Purification and characterization. *J Biol Chem.* 261:8754-8760.
- 10) Binkert C, Landwehr J, Mary JL, Schwander J, Heinrich G. 1989 Cloning, sequence analysis and expression of a cDNA encoding a novel insulin-like growth factor binding protein (IGFBP-2). *EMBO J.* 8:2497-2502.
- 11) Furlanetto RW, Underwood LE, Van Wyk JJ, D'Ercole AJ. 1977 Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay. *J. Clin Invest.* 60:648-657.
- 12) Daughaday WH, Kapadia M, Mariz I. 1987 Serum somatomedin binding proteins: physiologic significance and interference in radioligand assay. *J Lab Clin Med.* 109:355-363.
- 13) Breier BH, Gallaher BW, Gluckman PD. 1991 Radioimmunoassay for insulin-like growth factor-I: solutions to some potential problems and pitfalls. *J Endocrinol.* 128:347-357.
- 14) Daughaday WH, Mariz IK, Blethen SL. 1980 Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites: a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol-extracted serum. *J Clin Endocrinol Metab.* 51:781-788.
- 15) Blum WF, Gallaher B, Ranke MB. 1992 An IGFBP-blocked IGF-I RIA that measures what it pretends to measure: IGF-I. 74th Annual Meeting of the American Endocrine Society. 293.
- 16) Rosenfeld RG, Wilson DM, Lee PDK, Hintz RL. 1986 Insulin-like growth factors I and II in evaluation of growth retardation. *J Pediatr.* 109:428-433.
- 17) Clemmons DR, Van-Wyk JJ, Ridgway EC, Kliman B, Kjellberg RN, Underwood LE. 1979 Evaluation of acromegaly by radioimmunoassay of somatomedin-C. *N. Engl J Med.* 301:1138-1142

- 18) Zapf J, Walter H, Froesch ER. 1981 Radioimmunological determination of insulin-like growth factors I and II in normal subjects and in patients with growth disorders and extrapancreatic tumor hypoglycemia. J Clin Invest. 68:1321-1330.
- 19) Blum WF. 1992 Insulin-like growth factors and their binding proteins. In: Ranke MB, ed. Functional Endocrinologic Diagnostics in Children and Adolescence. Mannheim: J + J Verlag; 102-117.
- 20) Rieu M, Girard F, Bricaire H, Binoux M. 1982 The importance of insulin-like growth factor (somatomedin) measurements in the diagnosis and surveillance of acromegaly. J Clin Endocrinol Metab. 55:147-153.
- 21) Blum WF, Ranke MB, Bierich JR. 1988 A specific radioimmunoassay for insulin-like growth factor II: the interference of IGF binding proteins can be blocked by excess IGF-I. Acta Endocrinol (Copenh).118:374-380.
- 22) Wilson EM, Oh Y, Rosenfeld RG (1997) Generation and characterization of an IGFBP-7 antibody: Identification of 31 kD IGFBP-7 in human biological fluids and Hs578T human breast cancer conditioned media. J Clin Endocrinol Metab Vol 82, 4:1301-1303
- 23) Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2000) Significance of basal IGF-I, IGFBP-3 and IGFBP-2 measurements in the diagnostics of short stature in children. Horm Res 54:60-68
- 24) Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2001) Relevance of IGF-I, IGFBP-3 and IGFBP-2 measurements during GH treatment of GH-deficient and non-GH-deficient children and adolescents. Horm Res 55:155-124

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SUMMARY – MEDIAGNOST IGF-I ELISA E25

Reconstitution / Dilution of Reagents		
Standards A-E	Reconstitution in Sample Buffer PP	1 ml
Control Serum KS1	Reconstitution in Sample Buffer PP	500 µl
Control Serum KS2	Reconstitution in Sample Buffer PP	500 µl
Wash Buffer WP	dilute in A. dest. (eg. total volume of 50 ml in a graduated flask and fill up to 1000 ml)	1:20
Sample + Control Sera KS1 and KS2: dilute 1:100 in Sample Buffer PP, mix immediately, incubate at least for 15 min, max. 2h. Use 50 µl for each well in the assay.		
Before conducting the assay equilibrate all reagents to room temperature.		

Assay Procedure for Double Determinations:

Pipette	Reagent	Position
50 µl	Antibody Conjugate AK	in <u>all</u> wells used
50 µl	Sample Buffer PP (blank)	A1 and A2
50 µl	Standard A (0.5 ng/ml)	B1 and B2
50 µl	Standard B (2.5 ng/ml)	C1 and C2
50 µl	Standard C (6 ng/ml)	D1 and D2
50 µl	Standard D (12 ng/ml)	E1 and E2
50 µl	Standard E (18 ng/ml)	F1 and F2
50 µl	Control Serum KS1	G1 and G2
50 µl	Control Serum KS2	H1 and H2
50 µl	Samples	following wells
Cover the wells with the sealing tape.		

Incubation: 1 h at RT, 350 rpm

5x 300 µl	Aspirate the contents of the wells and wash 5x with 300 µl Wash Buffer WP	each well
100 µl	Enzyme Conjugate EK	each well

Incubation: 30 min at RT, 350 rpm

5x 300 µl	Aspirate the contents of the wells and wash 5x with 300 µl Wash Buffer WP	each well
100 µl	Substrate S	each well

Incubation: 30min in the dark RT

100 µl	Stop Solution SL	each well
Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength.		

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REF E25

International Test Description

CAL A-E	A -E	Rec in 1 ml BUF PP	
Control	KS1&KS2	Rec in 500 µl BUF PP	
WASHBUF 20x	WP		1:20 DILU A. dest.

SPE + Control 1:100 DILU BUF PP ↔ ⌚ 0.25 h
°C 20-25 °C

50 µl	Ab	A1 - End
50 µl	BUF PP	A1/2
50 µl	CAL A Std A (0.5 ng/ml)	B1/2
50 µl	CAL B Std B (2.5 ng/ml)	C1/2
50 µl	CAL C Std C (6 ng/ml)	D1/2
50 µl	CAL D Std D (12 ng/ml)	E1/2
50 µl	CAL E Std E (18 ng/ml)	F1/2
50 µl	CONTROL KS 1	G1/2
50 µl	CONTROL KS 2	H1/2
50 µl	SPE	
TAPE		

⌚ 1 h °C 20-25 ↔ 350 rpm

5x 300 µl	5x WASHBUF WP
100 µl	CONJ
TAPE	

⌚ 0.5 h °C 20-25 ↔ 350 rpm

5 x 300 µl	5x WASHBUF WP
100 µl	SUBST TMB S

⌚ 30 min °C 20-25 ☀

100 µl	H ₂ SO ₄ SL
MEASURE	