ampli-sRANKL

(EN) ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF FREE HUMAN sRANKL IN EDTA PLASMA, HEPARIN PLASMA, SERUM OR CELL CULTURE SUPERNATANTS
CAT. NO. BI-20452 12 X 8 TESTS

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Additional information on our products is available on our website.

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www.bmgrp.com

1) INTRODUCTION

SRANKL, soluble receptor activator of nuclear factor (NF)- κ B ligand (also: Osteoprotegerin ligand, OPGL) a member of the tumor necrosis factor (TNF) family, is the main stimulatory factor for the formation of mature osteoclasts and is essential for their survival. RANKL is produced by osteoblastic lineage cells and activated T lymphocytes. It activates its specific receptor RANK that is located on osteoclasts and dendritic cells. The effects are counteracted by OPG (compare BI-20402 OPG ELISA), which is secreted by various tissues and acts as an endogenous soluble receptor antagonist. As the serum concentration of sRANKL is usually quite low in normal samples an additional enhancement system is used in this ELISA.

Indications:

- Postmenopausal and senile Osteoporosis
- Glucocorticoid-induced Osteoporosis
- Diseases with locally increased resorption activity
- Arthritis
- Oncology

2) CONTENTS OF THE KIT

CONT	KIT COMPONENTS	QUANTITY
PLATE	Human recombinant OPG pre-coated microtiter strips in strip holder	12 x 8 tests
WASHBUF	Wash buffer concentrate 20x, natural cap	1 x 50 ml
AB	Goat polyclonal biotinylated anti sRANKL antibody, green cap, ready to use	1 x 11 ml
STD	Standards, (0; 0.125; 0.25; 0.5; 1; 2 pmol/l), white caps, lyophilised	6 vials lyophilised
CTRL	Control, yellow cap, lyophilised, exact concentration after reconstitution see label	1 vial lyophilised
CONJ	Conjugate, (streptavidin-alkaline phosphatase), amber cap, ready to use	1 x 22 ml
AMPLIFIER A	Inorganic salts and buffered enzyme solution with tetrazolium violet, blue cap, in dropper bottle, ready to use	1 x 13 ml
AMPLIFIER B	Stabilised NADPH solution, red cap, in dropper bottle, ready to use	1 x 13 ml
STOP	Stop solution, white cap, ready to use	1 x 7 ml

3) ADDITIONAL MATERIAL ADDED TO THE KIT

- 2 self-adhesive plastic films
- QC protocol
- Protocol sheet
- Instruction manual for use

4) MATERIAL AND EQUIPMENT REQUIRED BUT NOT SUPPLIED

- Precision pipettes calibrated to deliver 50 μ l, 100 μ l, 200 μ l, 300 μ l and disposable tips
- Distilled or deionised water
- Plate washer is recommended for washing, alternative multichannel pipette or manifold dispenser
- Refrigerator with 4°C (2-8°C)
- ELISA reader for absorbance at 490 nm
- Graph paper or software for calculation of results
- Incubator with horizontal plate shaker at 300-500 rpm

5) REAGENTS AND SAMPLE PREPARATION

All reagents of the kit are stable at 4°C (2-8°C) until the expiry date stated on the label of each reagent.

Sample preparation:

Collect venous blood samples by using standardized blood collection tubes for serum or plasma. We recommend performing plasma or serum separation by centrifugation as soon as possible (e.g. 20 min at 2000 x g, preferably at 4°C (2-8°C)). If this is not possible store the samples at 4°C (2-8°C) prior to centrifugation (up to one day). The acquired plasma or serum samples should be measured as soon as possible. For longer storage aliquot samples and store at -25°C or lower. All samples should undergo only 4 freeze-thaw cycles. Lipemic or haemolysed samples may give erroneous results. Samples should be mixed well before assaying. For cell culture supernatants the normal standard curve can be used.

Samples with results above the highest standard (2 pmol/l) should be diluted 1:10 with diluted WASHBUF and retested.

For further information on sample stability please visit our website www.bmgrp.com/products/assay_characteristics/sample_stability or contact our customer service by e-mail export@bmgrp.at or by phone +43/ 1/ 29107-45.

Reconstitute / Handling:


STD (Standards) and CTRL (Control): Each in 700 µl deionised or distilled water at room temperature (18-26°C) for 15 min. Reconstituted STD and CTRL are stable at -25°C or lower until expiry date on label, avoid freeze-thaw cycles.

WASHBUF (Wash buffer): Dilute the concentrate 1:20 (e.g. 50 ml WASHBUF + 950 ml distilled water). Crystals in the buffer concentrate will dissolve at room temperature. The diluted buffer is stable at 4°C (2-8°C) until expiry date stated on label. Use only diluted WASHBUF (Wash buffer) for the assay performance.


AMPLIFIER A/B: The amplifiers are supplied in dropper bottles ready to use. Contamination of the Amplifier with ubiquitous Alkaline Phosphatase (tears, sweat, etc.) may diminish the quality of results. Therefore we do not recommend to remove the dropper stoppings!

Recommendation: Squeeze out only the needed amount of amplifier A or B from dropper bottles, each into clean vials with a calibration (1 drop is approximately 50 µl) and pipette 100 µl per well. Do not put liquid back into dropper bottles.

6A) PRINCIPLE OF THE ASSAY:

6B) PRINCIPLE OF THE AMPLIFICATION

7) ASSAY PROTOCOL

All reagents and samples must be at room temperature (18-26°C) before use in the assay.

Mark position for BLANK/STD/SAMPLE/CTRL (Blank/Standard/Sample/Control) on the protocol sheet.


Take microtiter strips out of the aluminium bag, take a minimum of one well as Blank. Store unused strips with desiccant at 4°C (2-8°C) in the aluminium bag. Strips are stable until expiry date stated on the label.

- Pipette 100 µl STD/SAMPLE/CTRL (Standard/Sample/Control) in duplicate into respective wells, except blank.
- Add 100 µl AB (biotinylated anti sRANKL antibody, green cap) into each well, except blank, swirl gently.
- Cover tightly and incubate over night (18-24h) at room temperature (18-26°C) while on an automatic plate shaker (300-500 rpm).**
- Aspirate and wash wells 5x with 300 µl diluted WASHBUF (Wash buffer), remove remaining WASHBUF by hitting plate against paper towel after the last wash.
- Add 200 µl CONJ (Conjugate, amber cap) into each well, except blank.
- Cover tightly and incubate for 1 hour at room temperature (18-26°C).**
- Aspirate and wash wells 5x with 300 µl diluted WASHBUF (Wash buffer), remove remaining WASHBUF by hitting plate against paper towel after the last wash.
- Add 100 µl of Amplifier A (blue cap) into each well.
- Add 100 µl of Amplifier B (red cap) into each well.
- 10. Incubate for 45 +/-10 min at room temperature (18-26°C).**
The optimal OD for the highest standard in the assay should recover between 2.0 - 2.8 with a blank of <0.2. The suggested incubation time above was established during the development of the assay. Each laboratory should establish its own optimal incubation time based on their conditions (i.e. temperature, time, etc). After 30 minutes, the assay should be read periodically (at 5 minutes intervals) prior to adding stop solution to ensure that adequate colour development (>2.0 OD) is achieved.
- Add 50 µl STOP (Stop solution, white cap) into each well.
- Measure absorbance immediately at 490 nm with reference 630 nm, if available.

8) CALCULATION OF RESULTS

Read the optical density (OD) of all wells on a plate reader. Subtract the blank OD from the values of STD, CTRL and sample. Construct the standard curve from the OD values of the STD. Use commercially available software or graph paper. Obtain sample concentration from this standard curve. The assay was evaluated with 4PL algorithm. Different curve fitting methods need to be evaluated by the user.

Example STD-curve:

The quality control protocol supplied with the kit shows the results of the final release QC for each kit. Data for optical density obtained by customers may differ due to various influences and/or due to the normal decrease of signal intensity during shelf life. However, this does not affect validity of results as long as an OD of 1.00 or more is obtained for the standard with the highest concentration and the value of the CTRL is in range (target range see label).

9) ASSAY CHARACTERISTICS

Normal range:	Median: females = 0.37 pmol/l (n = 635); male = 0.46 pmol/l (n = 394) Normal values and precision data (chapter 10) were obtained with human serum samples. It is recommended to establish the normal range for each laboratory.
Standard range:	0; 0.125; 0.25; 0.5; 1; 2 pmol/l
Conversion factor pg/ml to pmol/l:	1 pg/ml = 0.05 pmol/l (MW: 20 kD)
Sample volume:	100 µl human EDTA plasma, Heparin plasma, serum or cell culture supernatant
Detection Limit:	(0 pmol/l + 3 SD): 0.02 pmol/l
Incubation time:	Overnight / 1 h / 30 min

For further information on assay characteristics please visit our website www.bmgrp.com/products/assay_characteristics or contact our customer service by e-mail export@bmgrp.at or by phone +43/ 1/ 29107-45.

10) PRECISION

Intra-Assay: 2 samples of known concentrations were tested 16 times in 1 assay

Inter-Assay: 2 samples of known concentrations were tested in 10 assays

Intra-Assay (n=16)		
Mean (pmol/l)	0.45	0.68
SD	0.04	0.06
CV%	8	9

Inter-Assay (n=10)		
Mean (pmol/l)	1.60	1.42
SD	0.09	0.05
CV%	6	3

11) TECHNICAL HINTS

- Do not mix or substitute reagents with those from other lots or sources.
- Do not mix stoppers and caps from different reagents or use reagents between lots.
- Do not use reagents beyond expiration date.
- Protect reagents from direct sunlight.
- Amplifier solutions must be colourless.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.

12) PRECAUTIONS

All test components of human source were tested with 3rd generation tests against HIV-Ab and HBsAg; and were found negative. Nevertheless, they should be handled and disposed as if they were infectious.

All liquid reagents contain 0.01% Proclin 300 as preservative. Avoid contact with skin and mucous membrane.

Proclin 300 is not toxic in concentrations used in this kit. It may cause allergic skin reactions – avoid contact with skin or eyes.

- Do not pipette by mouth.
- Do not eat, drink, smoke or apply cosmetics where reagents are used.
- Wear gloves, glasses and lab jacket while performing this assay.
- Sulfuric acid is irritating to eyes and skin. Avoid contact with skin and mucous. Irritations are possible - Flush with water if contact occurs!!

13) LITERATURE

- Bone loss in relation to serum levels of osteoprotegerin and nuclear factor-kappaB ligand: the Tromsø Study. L Jorgensen, A Vik, N Emaus, J Brox, J-B Hansen, E Mathiesen, and P Vestergaard. Osteoporos Int, 2010; 21(6): 931-938.
- Alendronate reduces osteoclast precursors in osteoporosis. P D'Amelio, A Grimaldi, MA Cristofaro, M Ravazzoli, PA Molinatti, GP Pescarmona, and GC Isaia. Osteoporos Int, 2009; doi: 10.1007/s00198-009-1129-1.
- The Effect of Zoledronic Acid on Serum Dickkopf-1, Osteoprotegerin, and RANKL in Patients with Paget's Disease of Bone. SA Polyzos, AD Anastasilakis, Z Efsthathiadou, M Kita, I Litsas, A Avramidis, G Arsos, E Moraidis, S Gerou, V Pavlidou, A Papatheodorou, and E Terpos. Horm Metab Res, 2009; 41(11): 846-850.
- ADA-deficient SCID is associated with a specific microenvironment and bone phenotype characterized by RANKL/OPG imbalance and osteoblast. AV Sauer, E Mrak, R J Hernandez, E Zacchi, F Cavani, M Casiraghi, E Grunebaum, CM Roifman, MC. Cervi, A Ambrosi, F Carlucci, M G Roncarolo, A Villa, A Rubinacci, and A Aiuti. Blood, 2009; 114, 3216-3226.
- Bone markers predict cardiovascular events in chronic kidney disease. A Fahrleitner-Pammer, J Herberth, SR Browning, B Obermayer-Pietsch, G Wimsberger, H Holzer, H Dobnig, and HH Malluche. J Bone Miner Res, 2008; 23(11): 1850-8.
- Serum osteoprotegerin concentrations are decreased in women with the polycystic ovary syndrome. H F Escobar-Morreale, JI Botella-Carretero, MÁ Martínez-García, M Luque-Ramírez, F Álvarez-Blasco, and JL San Millán. Eur J Endocrinol, 2008; 159: 225 - 232.

- Constitutional Thinness: Unusual Human Phenotype of Low Bone Quality. B Galusca, M Zouch, N Germain, C Bossu, D Frere, F Lang, MH Lafage-Proust, T Thomas, L Vico, and B Estour. *J Clin Endocrinol Metab.* 2008; 93(1):110-117.
- Serum levels of receptor activator of nuclear factor kappaB ligand (RANKL) in healthy women and men. K Kersch-Schindl, J Wendlova, S Kudlacek, A Gleiss, W Woloszczuk, and P Pietschmann. *Exp Clin Endocrinol Diabetes*, 2008; 116(8): 491-5.

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ASSAY PROTOCOL AND CHECKLIST

PREPARATION OF REAGENTS:

- ☐ Bring all reagents to room temperature (18-26°C).
- ☐ Prepare reagents and samples as instructed.
- ☐ Bring unused and prepared components to the storage temperature mentioned in the package insert.
- ☐ Take microtiter strips out of the alu bag and mark positions on the protocol sheet.

TEST PROCEDURE:

- ☐ Step 1) Pipette 100 µl STD/SAMPLE/CTRL (standard/sample/control) into each well, except blank.
- ☐ Step 2) Add 100 µl AB (biotinylated anti sRANKL ab, green cap) into each well except blank, swirl gently.
- ☐ **Step 3) Cover tightly and incubate over night (18-24h) at room temperature (18-26°C) while shaking on an automatic plate shaker (300-500 rpm).**
- ☐ Step 4) Aspirate and wash wells with 300 µl WASHBUF (Wash buffer) five times. Remove remaining buffer by hitting plate against paper towel.
- ☐ Step 5) Add 200 µl CONJ (Conjugate, amber cap) into each well except blank.
- ☐ **Step 6) Cover tightly and incubate for 1 hour at room temperature (18-26°C).**
- ☐ Step 7) Aspirate and wash wells with 300 µl WASHBUF (Wash buffer) five times. Remove remaining buffer by hitting plate against paper towel.
- ☐ Step 8) Add 100 µl Amplifier A (blue cap) into each well.
- ☐ Step 9) Add 100 µl Amplifier B (red cap) into each well.
- ☐ **Step 10) Incubate for 45 +/-10 min at room temperature (18-26°C).**

The optimal OD for the highest standard in the assay should recover between 2.0 - 2.8 with a blank of <0.2. The suggested incubation time above was established during the development of the assay. Each laboratory should establish its own optimal incubation time based on their conditions (i.e. temperature, time, etc). After 30 minutes, the assay should be read periodically (at 5 minute intervals) prior to adding stop solution to ensure that adequate color development (>2.0 OD) is achieved.
- ☐ Step 11) Add 50 µl STOP (Stop solution, white cap) into each well.
- ☐ Step 12) Measure absorbance immediately at 490 nm with reference 630 nm, if available.