

REF BIO82
96 Well Plate



ARGUTUS MEDICAL

Serum Collagen IV EIA

Enzyme Immunoassay

Instructions for Use

FOR RESEARCH USE ONLY
Not for use in Diagnostic Procedures

Manufactured by: Daiichi Fine Chemical Co., Ltd. Japan

TABLE OF CONTENTS

INTENDED USE	3
BACKGROUND	3
ASSAY PRINCIPLE	3
COMPONENTS	4
PRECAUTIONS	5
STABILITY AND STORAGE	6
ADDITIONAL MATERIALS REQUIRED	6
PREPARATION OF REAGENTS	7
SAMPLE HANDLING AND STORAGE	7
ASSAY PROCEDURE	8
CALCULATION OF RESULTS	8
LIMITATIONS OF USE	8
PERFORMANCE CHARACTERISTICS	9
EXAMPLE OF CALIBRATION CURVE	10
WARRANTY	10
OTHER ARGUTUS MEDICAL SMARTASSAYS	11
REFERENCES	11

INTENDED USE

The Argutus Medical Serum Collagen IV EIA provides a method for the quantitative determination of collagen IV in human serum. Please contact Argutus Medical for further information regarding the assay of collagen IV in other tissue fluids. The Argutus Medical Serum Collagen IV EIA is for research use only. Not for use in diagnostic procedures.

BACKGROUND

Chronic liver disease comprises a number of progressive disorders which culminate in liver cirrhosis and which are characterized by excessive deposition of collagen.

Although various types of collagen (type I, III, IV, V and VI) increase in the liver with the progression of fibrosis, type IV collagen, a constituent of the basement membrane, is particularly noteworthy for the following reasons: its serum level correlates with hepatic levels of collagen IV¹, serum levels of collagen IV fall in response to effective therapy¹ and it is the earliest type of collagen to be synthesized during experimental liver injury^{2, 3}. Serum collagen IV levels are elevated in a variety of liver diseases⁴⁻⁶, in particular, serum collagen levels have been found to be predictive of therapy response in Hepatitis C infection⁷, and to be sensitive indicators of therapy response in abstaining alcoholics¹.

ASSAY PRINCIPLE

The Argutus Medical Collagen IV EIA is designed for the assay of serum collagen IV. It is a solid phase one-step sandwich EIA. Collagen IV in the sample is bound simultaneously by a solid-phase monoclonal antibody and a monoclonal antibody-enzyme conjugate, each directed at different antigenic sites. This results in the collagen IV molecule being sandwiched between the solid phase and enzyme-labelled antibodies. After removing unbound enzyme-labelled antibody and sample, the plate is incubated with enzyme substrate. The resultant colour development is directly proportional to the amount of collagen IV in the sample.

COMPONENTS

1. Antibody coated Microassay plate: 12x8 well strips coated with IgG directed against human collagen IV.

PLA

READY TO USE
2. Collagen IV Calibrator: Purified collagen IV in phosphate buffer (pH 7.0) with bovine serum albumin. 1000 µg/L Stock Solution (1mL). Contains 0.015% Geneticin as preservative.

CAL

STOCK
3. Dilution Buffer: Phosphate buffer (pH 7.0) containing bovine serum albumin and horse serum (5mL). Contains 30 mg/L Proclin 300 preservative.

BUF	DIL	1X
-----	-----	----

READY TO USE
4. Conjugate: Anti-collagen IV mouse Fab' conjugated to horseradish peroxidase (20 mL). Contains 30 mg/L Proclin 300 as preservative.

CONJ	1X
------	----

READY TO USE
5. Wash Concentrate: 10x Conc. Phosphate buffer with Tween 20 (PBT), (2 bottles of 50 mL). Contain 30 mg/L Proclin 300 as preservative.

BUF	WASH	10X
-----	------	-----

CONCENTRATE
6. Substrate: Stabilised liquid TMB solution (15 mL).

SUBS	TMB
------	-----

READY TO USE
7. Stop Solution: 1M Sulphuric Acid (15 mL)

SOLN	STP
------	-----

READY TO USE
8. Plate Seal: 1 sheet
9. Instructions for use

INS

10. Uncoated microassay plate

PRECAUTIONS

SAFETY

- The Argutus Medical Urinary Collagen IV EIA kit is intended for use by qualified laboratory staff only.
- The kit contains material of human origin that has been tested and found to be negative for Hepatitis B surface antigen, Hepatitis C and HIV antibodies. However, since no test can provide complete assurance, treat all materials as potentially infectious.
- The Stop Solution contains sulphuric acid which is corrosive. Avoid contact with the skin and eyes. If contact occurs, rinse off immediately with water and seek medical advice.
- The Substrate contains TMB that may irritate the skin and mucous membranes. Any substrate which comes in contact with the skin should be rinsed off with water.
- Dispose of all clinical specimens, infected or potentially infected material in accordance with good laboratory practice. All such materials should be handled and disposed of as though potentially infectious.
- Residues of chemicals and kit components are generally considered as hazardous waste. All such materials should be disposed of in accordance with established safety procedures.
- Wear protective clothing, disposable latex gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
- Do not pipette materials by the mouth and never eat or drink at the laboratory workbench.

PROCEDURAL

- Do not use kit or individual reagents past their expiry date.
- Do not mix or substitute reagents from different kit lot numbers.
- Deviation from the protocol provided may cause erroneous results.
- Performing the assay outside the time and temperature ranges specified may produce invalid results. Assays not falling within the established time and temperature ranges must be repeated.
- Reagent delivery should be aimed at midpoint of the side of the wells, taking care not to scratch the side with the pipette tip.
- Do not allow the wells to dry at any stage during the assay procedure.
- Care must be taken not to contaminate components and always use fresh pipette tips for each sample and component.
- Do not use reagents that are cloudy or that have precipitated out of solution.
- Ensure Wash Concentrate is mixed thoroughly and no crystals remain before reconstitution.
- High quality distilled or deionised water is required for the Wash Solution. The use of poor quality or contaminated water may lead to background colour in the assay.
- Allow all reagents to come to room temperature (20-27°C) and mix well prior to use.

- Avoid leaving reagents in direct sunlight and/or above 2-8°C for extended periods.
- Always use clean, preferably disposable, glassware for all reagent preparation.
- Ensure that the bottom surface of the plate is clean and dry before reading.
- Before commencing the assay, an identification and distribution plan should be established.

STABILITY AND STORAGE

1. All kit reagents should be stored at 2-8°C and are stable as supplied until the expiry date shown.
2. Prepared Wash Solution (PBT) is stable for up to one month at 2-8°C.
3. Prepared Calibrator solutions should not be stored.
4. Plate assay wells should be stored in sealed bags with dessicant at 2-8°C until required for use. Return unused wells to the storage bag together with dessicant.

ADDITIONAL MATERIALS REQUIRED

1. 20 µL and 150 µL micropipettes and a 100-150 µL multichannel pipette
2. Microassay strip washing system
3. ELISA plate reader capable of measuring at 450nm with reference at 630nm if available
4. 1 L beaker
5. Timer
6. Liquid trough
7. Deionised/Distilled water
8. Graduated cylinder (500 mL)

PREPARATION OF REAGENTS

WASH SOLUTION (PBT)

Perform a 1/10 dilution of Wash Concentrate adding, for example, 10 mL Wash Concentrate to 90 mL deionised water as required. Prepare only the volume of Wash Solution required for the assay. Each row of assay wells requires 10 mL of Wash Solution.

Ensure salt crystals are dissolved prior to dilution.

Gentle warming of Wash Concentrate at 37°C for 30 minutes will aid dissolution of salt crystals.

CALIBRATORS

Using labelled tubes prepare Calibrators as follows:

Collagen IV Concentration (µg/L)	Calibrator Volume (µL)	Dilution Buffer (µL)
1000 (A)	150 (A) (Stock Solution)	-
500 (B)	150 (A)	150
250 (C)	150 (B)	150
125 (D)	150 (C)	150
62.5 (E)	150 (D)	150
31.2 (F)	150 (E)	150
15.6 (G)	150 (F)	150
0 (H)	-	150

Calibrators should be prepared immediately before use. Do not store. The diluted calibrators are stable for at least 6 hours at 2-8° C.

SAMPLE HANDLING AND STORAGE

Samples can be stored at 2-8°C for one week. Samples may be stored at –20°C for 12 months. Repeated freeze-thawing of samples should be avoided.

ASSAY PROCEDURE

NOTE: All reagents should be allowed to reach room temperature prior to commencement of assay.

1. MIXING OF CALIBRATOR/SAMPLE

- 1.1 Prepare Wash Solution and Calibrators as described in "Preparation of Reagents".
- 1.2 Add Calibrators (**H-A: 0 – 1000 µg/L**) and samples (**20 µL/well**), in duplicate, to the uncoated Vinyl Microassay plate.
- 1.3 Add **150 µL** Conjugate to each well.

2. IMMUNOREACTION

- 2.1 Place required number of anti-collagen IV coated Microassay wells in the assay plate (16 for the Calibrators plus two each for the samples).
- 2.2 Transfer **100µL** of the mixtures from above into the equivalent wells in the anti-Collagen IV coated Microassay wells.
- 2.3 Cover the Microassay plate with the lid and incubate at room temperature (20-27°C) for exactly **30 minutes**.
- 2.4 Remove the plate lid and wash each strip three times (**350 µL/well**) with Wash Solution. When complete, firmly tap the plate against a paper towel to ensure complete removal of wash fluid from wells.

3. COLOUR DEVELOPMENT

- 3.1 Add **100 µL** Substrate/well using a multichannel pipette and incubate at room temperature (20-27°C) for exactly 30 minutes.

4. STOP

- 4.1 Add **100 µL** Stop Solution/well using a multi channel pipette. Ensure complete mixing of Substrate and Stop Solution.
- 4.2 Read **immediately** at 450nm using 630nm as reference (if available).

CALCULATION OF RESULTS

1. Calculate the mean absorbance for each calibrator and sample.
2. Plot a Calibration curve of $A_{450/630nm}$ versus collagen IV (µg/L) on a log-log scale.
3. Read the collagen IV (µg/L) indicated by the mean absorbances of the samples from the calibration curve.
4. If the samples have been diluted, multiply the calculated [collagen IV] by the appropriate dilution factor in order to obtain the actual [collagen IV].
5. Concentrations of samples with readings outside the standard curve are invalid and must be repeated with a higher dilution factor. It is not acceptable to extrapolate data.

PERFORMANCE CHARACTERISTICS

NORMAL RANGE

Based on healthy Japanese volunteers, the reference normal range for Collagen IV is: $99 \pm 23 \mu\text{g/L}$. Mean \pm 1S.D. (N = 180).

LIMIT OF DETECTION

The detection limit of the Argutus Medical Collagen IV EIA is $15.6 \mu\text{g/L}$.

MEASURING RANGE

The calibration curve range covers the range $15.6\text{--}1000 \mu\text{g/L}$. This range may be extended by increasing sample dilution.

SPECIFICITY

The Argutus Medical Serum Collagen EIA IV is highly specific for the detection of collagen IV. Cross reactivity is less than 2% with Collagen II and less than 0.5% with other forms of collagen.

SENSITIVITY

When reading from the standard curve the $A_{450\text{nm}}$ value of the $1000 \mu\text{g/L}$ standard should be >0.6 OD.

INTERFERENCE

No significant interference has been observed in this assay with lipaemic, haemolytic or icteric samples.

- Lipaemia: Less than 10% interference up to 1200 Formazine turbidity units.
- Haemolysis: Less than 10% interference up to 3 g/L haemoglobin.
- Icterus: Less than 10% interference up to 0.2 g/L bilirubin.

DILUTION - RECOVERY

Dilution of samples containing high levels of collagen IV gave the following results:

Sample	Dilution								
	1 / 2			1 / 4			1 / 8		
	Expected $\mu\text{g/L}$	Obtained $\mu\text{g/L}$	Recovery %	Expected $\mu\text{g/L}$	Obtained $\mu\text{g/L}$	Recovery %	Expected $\mu\text{g/L}$	Obtained $\mu\text{g/L}$	Recovery %
A	57	61	107	28	32	114	14	16	114
B	107	110	103	53	58	109	27	28	104
C	259	270	104	130	139	107	65	67	103

REPRODUCIBILITY

Intra-assay variation of the Argutus Medical Collagen IV EIA

Sample	\bar{x} [Collagen IV] $\mu\text{g/L}$	SD	%CV	N
Low	119	7.4	6.2	8
Medium	218	7.8	3.6	8
High	520	12	2.3	8

Inter-assay variation of the Argutus Medical Collagen IV EIA

Sample	\bar{X} [Collagen IV] $\mu\text{g/L}$	SD	%CV	N
Low	115	11	9.6	6
Medium	291	13	4.5	6
High	370	31	8.2	6

Inter-batch Variation of the Argutus Medical Collagen IV EIA calculated for three batches of kits

Sample	\bar{X} [Collagen IV] $\mu\text{g/L}$	SD	%CV	N
Low	115	5	4.3	3
Medium	270	6	2.2	3
High	386	16	4.1	3

EXAMPLE OF CALIBRATION CURVE




Figure 1: Typical Calibration curve obtained using the Argutus Medical Collagen IV EIA.
Plot of $A_{450/630\text{nm}}$ versus [Collagen IV] $\mu\text{g/L}$.
Assay range is 15.6 – 1000 $\mu\text{g/L}$.

WARRANTY

The performance data presented here were obtained using the procedure described. Any change or modification of the procedure, not recommended by Argutus Medical, may affect the results, in which case Argutus Medical disclaims all warranties, expressed, implied or statutory, including implied merchantability and fitness for use. In the case of such an event, Argutus Medical shall not be liable for damages, direct or consequential.

OTHER ARGUTUS MEDICAL SMARTASSAYS

Pancreatic Injury Testing Service

Catalogue No	Product Name	Description
TEST BBU	Trypsinogen Activation Peptide (TAP) EIA	TAP in human and mammalian urine and tissue

Animal Organ Damage Biomarkers

Catalogue No	Product Name	Description
BIO64RT	Rat Alpha GST EIA	α GST in rat serum, urine and tissue culture
BIO76YB1	Rat Yb1 GST EIA	GSTYb1 (μ GST) in rat urine
BIO89RPA1	RPA-1 EIA	Renal papillary antigen 1 in rat urine
BIO87CD	RPA-1 Antibody	Antibody to rat collecting duct
BIO88LH	RPA-2 Antibody	Antibody to rat loop of henle

Human Organ Damage Biomarkers

Catalogue No	Product Name	Description
BIO66NEPHA	NEPHKIT® Alpha GST EIA	α GST in human urine
BIO60HEPA	HEPKIT® Alpha GST EIA	α GST in human serum and plasma
BIO60HEPAS	High Sensitivity Alpha GST EIA	α GST in human serum and plasma
BIO85	PI GST EIA	π GST in human urine and plasma
BIO83	Urinary Collagen IV EIA	Collagen IV in human urine
BIO81DNA	OxyDNA test	Fluorescence method for the detection of oxidative DNA damage in cell suspensions

REFERENCES

- 1) **Tsutsumi, M. et al.** (1996). Serum biomarkers for hepatic fibrosis in alcoholic liver disease: which is the better biomarker, type III procollagen, type IV collagen, laminin, tissue inhibitor of metalloprotease or prolyl hydroxylase? *Alcoholism: Clinical and Experimental Research*. **20(9)**, 1512-1517.
- 2) **Gay, S. and Miller, E.J.** (1983). What is collagen, what is not. *Ultrastruct. Pathol.* **4**, 365-377.
- 3) **Dingleman, R.F. et al.** (1983). Collagen formation by the hepatocyte in primary monolayer culture and *in vivo*. *Science*, **219**, 1343-1345.
- 4) **Obata, K. et al.** (1989). One step sandwich enzyme immunoassay for human type IV collagen using monoclonal antibodies. *Clin. Chim. Acta*, **181**, 293- 304.
- 5) **Yokoya, Y. et al.** (1992). Concentration of serum laminin and type IV collagen in liver diseases assayed by a sandwich enzyme-immunoassay using monoclonal antibodies. *Clin. Chim. Acta*, **210**, 109- 118,
- 6) **Ueno, T. et al.** (1992). Significance of serum type-IV collagen levels in various liver diseases. *Scandinavian Journal of Gastroenterology* **27**, 513-520.
- 7) **Yabu, K. et al.** (1994). Serum collagen IV for the assessment of fibrosis and resistance to Interferon therapy in chronic hepatitis C. *Scandinavian Journal of Gastroenterology*. **29**, 474-479.



ARGUTUS MEDICAL

**Argutus Medical Ltd.
Unit 9 Trinity Technology & Enterprise Campus,
Pearse Street,
Dublin 2,
Ireland**

Tel: +353 1 670 8576

Fax: +353 1 670 8575

info@argutusmed.com

<http://www.argutusmed.com>

Manufactured by: Daiichi Fine Chemical Co, Ltd, Japan.

Colls-437-05 03/09