# Resistin

# ELISA

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## PACKAGE INSERT ENGLISH

### Mediagnost Resistin ELISA E50

- is suited for Resistin determination in **Serum** and **Plasma** samples
- is extremely **sensitive** (**12 pg/ml  $\cong$  1.2 pg per well**) and, thus allows measurements in cell culture media too and in specimens others than serum e.g. in Cerebrospinal fluid, Amnion fluid, Saliva, Urine, Breast milk
- is **fast**: incubation time a total of 4 hours
- Single Standards with **20, 100, 300, 600,1000 pg/ml** human Resistin are provided in the Kit
- Control Serum **KS** is of human serum
- is calibrated with **recombinant Resistin**
- Microtiter plates are separately breakapart, tests can be adapted to individual requirements

### INTENDED USE

Measurement of human Resistin in human Serum and Plasma Sample

### CLINICAL IMPLICATION

- Resistin is relevant e.g. in research of:
- Adiposity
- Insulin Resistance, Diabetes
- Arteriosclerosis
- Inflammation

### INTRODUCTION

Resistin, a cysteine-rich protein of 11.3 kDa (1), was firstly found in mice (2) and constitutes together with RELM $\alpha$ , RELM $\beta$  and RELM $\gamma$  the protein family of resistin-like molecules (RELM).

In humans, resistin and RELM $\beta$  (1) but no other proteins of the RELM family were found. The human form of resistin shows a homology of 53% to the murine protein (4). It has 11 cysteine-residues, is synthesized as a propeptide of 108 amino acids and secreted as a dimer, build by a disulfide bridge of cysteine residues (22). Beside this intermolecular disulfide bridge, 5 additional intramolecular ones exist (5,6).

Appearance of multi- and oligomer formation was proved by size exclusion chromatography. Thereby it was shown, that oligomer formation is SDS-insensitive but can be inhibited by  $\beta$ -mercaptoethanol and is therefore likely to be caused by disulfide bridges (1). Further on, the resistin structure seems to be dependent on its concentration, as circular dichroism analysis shows a concentration dependent shift of  $\alpha$ -helical to  $\beta$ -sheet structure (1).

Resistin expression was demonstrated in white adipose tissue (10), pituitary (11) and pancreatic islets (12) of mice as well as in brown adipose tissue of rats. In humans, resistin expression in adipocytes can be detected but only at a very low level. But in vitro, resistin expression of non-adipocytes in fatty tissue was shown (13). Human resistin gene is also expressed in pancreatic islets (12), pre-adipocytes (14) macrophages (15) and bone marrow (39). So, resistin is of relevance for inflammation processes as well as for lipid metabolism.

Most investigation refers to the mouse model. Here, the existence of trimeric and hexameric resistin in serum was demonstrated (7). In comparison to adiponectin biology it is highly probable that different resistin oligomers have different biologic function (8, 9).

In mice, a correlation between adiposity, insulin resistance and resistin expression was found empirically. In humans, respective study results are not clear – several studies show an association of resistin serum concentration and adiposity or insulin resistance (17, 25-31). But others failed in confirming these results (14, 16-24). Therefore, there is requirement for valid and reproducible determination of resistin serum concentration.

Relevance of resistin in other physiologic processes than energy metabolism was investigated by several different approaches. Experiments with endothelial cells gave interesting results. Here, resistin was shown to enhance expression of VCAM-1 and ICAM-1 (33, 34). By this way, resistin is potentially able to influence endothelial inflammation (35, 36) and, thereby atherosclerosis. These results were confirmed by experiments in mice, where endothelin-1 was shown to regulate resistin secretion (37, 38).

In recent research human resistin was shown to increase pre-adipocyte proliferation and lipolysis of mature adipocytes (38). By the way of modulating MAPK-signalling pathways resistin exerts crucial influence on energy metabolism.

Present research demonstrates, that Resistin exerts influence on a broad variety of physiological processes, however a clear and defined biological role of resistin remains still unexisting.

This ELISA-kit enables the user to determine the exact concentration of Resistin in human serum/plasma as well as other body fluids and thereby assists investigation of Resistin biology.

## REAGENTS PROVIDED

1)	MTP	<b>Microtiter plate</b> , ready for use, with 96 wells, dived up in 12 stripes à 8 wells (separately breakapart), coated with human Resistin antibody.
2)	CAL	<b>Standards A-E, lyophilised, contain recombinant Resistin. Standard values are between 0.02 - 1 ng/ml (20, 100, 300, 600 und 1000 pg/ml) Resistin and have to be reconstituted with 750 µl (each) Sample Buffer PP.</b> Attention: Please use only Sample Buffer PP for this dilution, because only this assures, that the Standards and the respective samples subsequently will incubate under identical conditions in the same special buffer!
3)	BUF	<b>Sample Buffer PP, 120 ml</b> , ready for use, please use for the reconstitution of the Standards A – E and for the sample and Control Serum KS dilution.
4)	BUF	<b>Dilution buffer VP, 25 ml</b> , ready for use, please use this for the <b>reconstitution of Control Serum KS</b> and for the <b>dilution of Antibody Conjugate AK and Enzyme Conjugate EK.</b>
5)	Control	<b>Control Serum KS</b> , lyophilised: Contains human Serum and has to be reconstituted with <b>100 µl Dilution buffer VP.</b> The Resistin target value concentration and the respective range are given on the vial label. The dilution of the <b>Control Serum KS</b> in <b>Sample Buffer PP</b> should be according to the dilution of the respected samples.
6)	Ab	<b>Antibody Conjugate AK, 120 µl, 100fold concentrated solution</b> , contains biotinylated anti-Resistin antibody, please dilute before use 1:100 in <b>Dilution buffer VP</b> : e.g., add 100 µl <b>Antibody Conjugate AK</b> to 10 ml <b>Dilution Buffer VP</b> , mix and use 100 µl/well of this dilution in the assay.
7)	CONJ	<b>Enzyme Conjugate EK, 120 µl, 100fold concentrated solution</b> , contains HRP (Horseradish peroxidase)-labelled Streptavidin, please dilute before use 1:100 in Dilution Buffer VP: e.g. add 100 µl Enzyme conjugate EK to 10 ml Dilution buffer VP, mix and use 100 µl/well of this dilution in the assay.
8)	WASHBUF 20x	<b>Washing Buffer (WP), 50 ml, 20fold concentrated solution. Washing Buffer (WP)</b> has to be diluted 1:20 with distilled or demineralised water before use (e.g. add the complete contents of the flask (50 ml) into a graduated flask and fill up with A. dest. to 1000 ml). Attention: After dilution the Washing Buffer is only 4 weeks stable at 2-8°C, dilute only according to requirements.
9)	SUBST	<b>Substrate (S), 12 ml</b> , ready for use, horseradish-peroxidase-(HRP)-substrate, stabilised H <sub>2</sub> O <sub>2</sub> Tetramethylbencidine.
10)	H <sub>2</sub> SO <sub>4</sub>	<b>Stopping Solution (SL), 12 ml</b> , ready for use, 0.2 M sulphuric acid, Caution acid!
11)		Sealing tape for covering of the microtiter plate, 3 x, adhesive.

## MATERIALS REQUIRED BUT NOT PROVIDED

Precision pipettes and multichannel pipettes with disposable plastic tips  
Distilled or deionized water for dilution of the Washing Buffer (WP)  
Microtiter plate shaker (350 rpm)  
Microtiter plate washer (recommended)  
Micro plate reader ("ELISA-Reader") with filter for 450 ≥590 nm.  
Polyethylene PE/Polypropylene PP tubes for dilution of samples

## WARNINGS AND PRECAUTIONS

### For In Vitro Diagnostic Use only. For professional use only.

Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.

Before use, all kit components should be brought **to room temperature at 20 - 25°C**. Precipitates in buffers should be dissolved before use by thorough mixing and warming. **Temperature will affect the absorbance** readings of the assay. However, values for the patient samples will not be affected.

Do not mix reagents of different lots. Do not use expired reagents.

Disposal of containers and unused contents should be done in accordance with federal and local regulatory requirements.

Do not use obvious damaged or microbial contaminated or spilled material.

Caution: This kit contains material of human and/or animal origin. Source human serum for the Control Serum provided in this kit was tested by FDA recommended methods and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV) antibodies. No known test methods can offer total assurance of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

### 2-Methyl-4-Isothiazolin-3-one

Following components contain < 0.01% **2-Methyl-4-isothiazolin-3-one** solution as preservative **A-E, AK, EK, VP, PP** < 0.01% 2-Methyl-4-isothiazolin-3-one Solution

R36/38 Irritating to eyes and skin

R43 Sensibilisation through skin contact possible

S26 In case of contact with eyes rinse immediately with plenty of water and seek medical advice

S28.1 After contact with skin wash immediately with plenty of water

### 5-chloro-2-methyl 2H isothiazol-3-one and 2-methyl-2H-Isothiazol-3-one

Following components contain < 0.01%(w/w) 5-chloro-2-methyl 2H isothiazol-3-one and 2-methyl-2H-isothiazol-3-one as preservative: **A-E, AK, EK, VP, PP, WP**

R36/38 Irritating to eyes and skin

R43 Sensibilisation through skin contact possible

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice S28.1

S28.1 After contact with skin, wash immediately with plenty of water

### Stop solution contains 0.2 M Sulfuric Acid (H<sub>2</sub>SO<sub>4</sub>)

R36/38 Irritating to eyes and skin

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

S28.1 After contact with skin, wash immediately with plenty of water

S36/37 Wear suitable protective clothing and gloves.

Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step. Use separate pipette tips for each sample, control and reagent to avoid cross contamination. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.

**TMB-Substrate (S) contains 3,3',5,5' Tetramethylbenzidine.** Store and incubate in the dark.

R20/21/R22 Harmful by inhalation, in contact with skin and if swallowed

R36/37/38 Irritating to eyes, respiratory system and skin

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

S28.1 After contact with skin, wash immediately with plenty of water

S36/37 Wear suitable protective clothing and gloves

### General first aid procedures:

Skin contact: Wash affected area thoroughly with water. Discard contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: If swallowed, wash out mouth thoroughly with water. Immediately see a physician.

Do not eat, drink or smoke in these areas.

Never pipette the materials with the mouth.

Spilled material must be wiped off immediately and should become disinfected. Clean contaminated areas and equipment with a suitable detergent.

## METHOD

The enzyme immunoassay for Resistin E50 is a so-called Sandwich-Assay. It utilizes a specific high affinity polyclonal rabbit antiserum coated on the wells of a microtiter plate. The Resistin in the samples binds quantitatively to the immobilized antiserum. In the following step, the biotinylated antiserum binds in turn to Resistin. After washing, Streptavidin-Peroxidase-Enzyme conjugate will be added, which will bind highly specific to the biotin of the antiserum and will catalyse in the closing substrate reaction the turn of the colour, quantitatively depending on the Resistin level of the samples.

## SPECIMEN

Serum as well plasma samples are suitable (significant deviation of Resistin levels in corresponding serum-, Heparin-, EDTA-, Citrate-plasma-Samples were not found). Haemolytic samples appear to show falsely high Resistin levels, using such samples should be checked out critically. Common cell culture medium, saliva, breast milk and urine were found to be suitable specimens too.

By means of the special sample buffer an external sample preparation prior to the assay is not required (see below).

The blood sample for serum preparation should be gained according to standardized venipuncture procedure. The samples should be stored without anticoagulation reagents. Haemolytic reactions have to be avoided. The blood has to be allowed to clot and after complete clotting, serum is separated by centrifugation.

## Storage of the samples

Storage at RT           max. 2 days

Storage at  $-20^{\circ}\text{C}$    max. 2 years

in tightly closable plastic tubes.

More than 3 freeze/thaw cycles are not possible.

## Sample Preparation

Samples have to be diluted in Sample Buffer (PP). The excellent linearity of this test system allows sample dilution of 1:5 to 1:400.

In most determinations (serum or plasma samples, and no extreme values expected) a dilution **from 1:10 to 1:50 with Sample Buffer PP** should be suitable. According to expected Resistin levels the dilution with PP can be higher or lower. Because the sample buffer **PP** is special composed for the correct determination of Resistin, the dilution should be **at least 1:5!**

Resistin concentrations may be completely different in body fluids of human origin other than serum or cell culture supernatants.

For clinical purposes we recommend a standard dilution of 1:21.

#### Suggestion for dilution protocol:

Pipette 300 µl Sample Buffer PP in PE-/PP-Tubes (application of a multi-stepper is recommended in larger series), add 15 µl Serum- or Plasma (dilution 1:21). After mixing use 2 x 100 µl of this dilution in the assay.

## **ASSAY PREPARATION AND TECHNICAL NOTES**

The assay has to be conducted strictly according the test protocol herein.

Reagents with different lot numbers cannot be mixed.

The microtiterplate and reagents are stable until the indicated expiry if stored unopened and protected from sunlight at 2 – 8°C.

The **shelf life** of the components after initial opening is limited to **4 weeks**, if stored appropriately.

**Bring all reagents to room temperature (20 - 25°C) before use.** Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming.

### **Incubation at room temperature means: 20-25°C**

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtiter plate shaker. We are recommending 350 rpm. Due to certain technical differences deviations may occur, in case the rotation frequency must become adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/or false values, excessive shaking may result in high optical densities and/or false values.

Proper washing is of basic importance for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided washing buffer diluted to usage concentration. Washing volume per washing cycle and well must be 300 µl at least.

The danger of handling with potentially infectious material must be taken into account.

When using an automatic microtiter plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non fuzzy absorbent tissue.

Manual washing is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamically swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non fuzzy absorbent tissue.

### **Standards and Control**

For the reconstitution of the lyophilised **Standards A - E Sample Buffer PP** has to be used.

The lyophilised **Control Serum KS** must be reconstituted with the **Dilution Buffer VP**. The dilution of the **Control Serum in Sample Buffer PP** should be according the dilution of the respected samples.

It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam!) with a Vortex mixer.

The reconstituted standard and controls can be stored for **4 weeks at -20°C**. Repeated freeze/thaw cycles have to be avoided.

### Antibody and Enzyme Conjugate

Use the Dilution Buffer **VP** for the dilution of the Antibody Conjugate **AK** and Enzyme Conjugate **EK** 100fold concentrates. The diluted solutions are only limited stable at 2-8°C and should be prepared daily fresh.

### Washing Buffer

The required volume of washing buffer is prepared by 1:20 dilution of the provided 20fold concentrate with deionised water. The diluted Washing Buffer is stable for 4 weeks at 2-8°C. It has to be at room temperature for usage!

### Microtiterplate

After initial opening, store the unused strips and microtiter wells together with the desiccant at 2-8°C in the clip-lock bag.

### Substrate Solution

The Substrate Solution (S), stabilised H<sub>2</sub>O<sub>2</sub>-Tetramethylbencidine, is photosensitive – store and incubate in the dark.

## ASSAY PROCEDURE

NOTES: All determinations (Standards, Control Serum and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended. When performing the assay, the Standards, Control Serum and the samples should be pipette as fast as possible (e.g., <15 minutes). To avoid distortions due to differences in incubation times, **Antibody-Conjugate AK** and the **Enzyme Conjugate EK** as well as the following **Substrate Solution S** should be added to the plate in the same order and in the same time interval as the samples. **Stop Solution SL** should be added to the plate in the same order as the Substrate Solution.

- 1) Add **100 µl Sample Buffer PP** in wells A1/A2 (blank) and
- 2) Pipette in positions B1/2 **100 µl of the Standard A** (0.02 ng/ml),  
pipette in positions C1/2 **100 µl of the Standard B** (0.1 ng/ml),  
pipette in positions D1/2 **100 µl of the Standard C** (0.3 ng/ml),  
pipette in positions E1/2 **100 µl of the Standard D** (0.6 ng/ml),  
pipette in positions F1/2 **100 µl of the Standard E** (1 ng/ml).
- 3) To control the correct accomplishment **100 µl** of the 1:21 (or in respective dilution rate of the sample) in **Sample Buffer PP** diluted **Control Serum KS** can be pipetted in positions G1/2.
- 4) Pipette **100 µl** each of the **diluted sample** (e.g. dilute 1:21 with Sample Buffer **PP**) in the rest of the wells, according to requirements.
- 5) Cover the wells with sealing tape and incubate the plate for **2 hours** at **room temperature** (shake at 350 rpm). After incubation aspirate the contents of the wells and wash the wells 5 times with **300 µl Washing buffer WP** / well.
- 6) Following the last washing step pipette **100 µl** of the 1:100 with **Dilution Buffer VP** diluted **Antibody Conjugate AK** in each well. Cover the wells with the sealing tape and incubate **1 hour** at room temperature (shake at 350 rpm).
- 7) After incubation wash the wells 5 times with **Washing Buffer WP** as described in step 5).

- 8) Following the last washing step, pipette **100 µl** of the 1:100 with Dilution Buffer VP diluted **Enzyme Conjugate EK** in each well. Cover the wells with the sealing tape and incubate for **30 minutes** at room temperature (shake at 350 rpm).
  - 9) After incubation wash the wells 5 times with Washing Buffer **WP** as described in the step 5.
  - 10) Pipette **100 µl** of the **TMB-substrate** solution **S** in each well.
  - 11) Incubate the plate for **30 minutes** in the dark at **room temperature**.
  - 12) Stop the reaction by adding **100 µl** of **Stopping Solution SL** to all wells.
- Measure the absorbance within **30 minutes** at **450 nm** (reference filter: **≥ 590 nm**)

## ESTABLISHING THE STANDARD CURVE

For the evaluation of the assay it is preconditioned that the absorbance values of the blank should be below 0.3, these of standard E should be above 0.8.


Samples, which yield higher absorbance values than Standard E are beyond the standard curve, for reliable determinations these samples should be tested anew with a higher dilution.

The standards provided contain the following concentrations of Resistin:

<b>Standard</b>	<b>A</b>	<b>B</b>	<b>C</b>	<b>D</b>	<b>E</b>
ng/ml	0.02	0.10	0.30	0.60	1.00
pg/ml	20	100	300	600	1000

- 1) Calculate the mean absorbance value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the blank from the mean absorbances of all other values.
- 3) Plot the standard concentrations on the x-axis versus the mean value of the absorbance of the standards on the y-axis.
- 4) Recommendation: Calculation of the standard curve should be done by using a computer program because the curve is in general (without respective transformation) not ideally described by linear regression. **A higher-grade polynomial, or four parametric logistic (4-PL) curve fit or non-linear regression** are usually suitable for the evaluation (as might be spline or point-to-point alignment in individual cases).
- 5) The **Resistin concentration** of the diluted sample or the diluted control serum KS in ng/ml (or µg/ml according the chosen unit for the standards) is calculated in this way, the Resistin concentration of the **undiluted sample** and of KS is calculated **by multiplication with the respective dilution factor**.

### E50 Standard Curve

**Fig. 1. Exemplary Standard Curve** with a polynomial 2<sup>nd</sup> degree as curve fit.

The exemplary shown standard curve in Fig.1 **cannot** be used for calculation of your test results. You have to establish a standard curve for each test you conduct!

Exemplary calculation of the Resistin concentration of a 1:21 diluted sample:

Measured extinction of your sample	0.85
Measured extinction of the blank	0.05

Your measurement program will calculate the Resistin concentration of the diluted sample automatically by using the difference of sample and blank for the calculation. You only have to determine the most suitable curve fit (here: polynomial 2<sup>nd</sup> degree).

In this exemplary case the following equation is solved by the program to calculate the Resistin concentration in the sample:

$$y = 0.0046207 + 3.3674 x - 1.5709 x^2$$

$$0.2686 = x$$

if the dilution factor (1:21) is taken into account the Resistin concentration of the undiluted sample is  $0.2686 \times 21 = 5,64 \text{ ng/mL} = 0.00564 \text{ } \mu\text{g/mL}$

## PERFORMANCE CHARACTERISTICS

### Standards

The standards are prepared from recombinant human Resistin (19.5 kDa, 2 x 92 amino acids, expressed in *E. coli*) in concentrations of 20, 100, 300, 600 and 1000 pg/ml (pico Gramm / ml, equal to 0.02 ng/ml-1 ng/ml).

### Sensitivity

The **analytical sensitivity** of the assay yields **0.012 ng/ml** (12 pg/ml; as 2x SD of zero standard in 15fold determination).

### Specificity

Commercially available sera from bovine, cat, chicken, dog, donkey, goat, guinea pig, horse, mouse, pig, rabbit, rat and sheep were diluted (1:10) and used as samples in this assay system and the signal intensity was measured. No cross reactivity was detected.

### Interference

Interference of physiological appearing substance with the Resistin measurement was investigated. Serum samples have been enriched with different concentrations of possibly interfering substances and the amount of Resistin was measured and compared with the Resistin concentration in the same sample without any enrichment. In table 1 the relative results are shown. None of the tested substances interfered significantly with Resistin measurement.

**Table 1:** Interference: Three serum samples were enriched with indicated amount of the potentially interfering substance and measured. Shown is % of Resistin of the native, non enriched serum sample

	<b>Triglyceride 100 mg/ml</b>	<b>Bilirubin 100 µg/ml</b>	<b>Haemolysate 1000 µg/ml</b>
Serum 1	101	93	94
Serum 2	115	99	99
Serum 3	104	103	147

**Table 2:** Effects of coagulation inhibitors were investigated by adding indicated amounts of inhibitors to PP enriched with 0.3 ng/ml Resistin. Relative amounts of Resistin measured in inhibitor containing samples in comparison to 0.3 ng/ml Resistin containing Sample Buffer (PP) are shown.

		<b>% of Resistin in PP</b>	
		<b>Mean (n=3)</b>	<b>SD</b>
3.8 g/l	Citrate	94	7.67
0.0068 mol/l	EDTA	93	4.96
30,000 IE/l	Heparin	96	4.89

## Reproducibility and Precision

The inter- and intra assay coefficients of variability are below than 6.8% and 5%, respectively. Exemplary determinations are shown in table 3 and table 4.

**Table 3:** Intra-Assay-Variation

	Number of determinations	Mean value (µg/ml)	Standard deviation (µg/ml)	VC (%)
Sample 1	16	5.87	0.138	2.35
Sample 2	16	12.19	0.377	3.10
Sample 3	6	14.36	0.668	4.66

**Table 4:** Inter-Assay-Variation (results of 11 determinations, each)

	Mean value (ng/ml)	Standard deviation (ng/ml)	VK (%)
Sample 1	2.70	0.16	5.94
Sample 2	4.20	0.28	6.77
Sample 3	5.80	0.28	4.79

## Recovery and Linearity

The Mediagnost Resistin ELISA E50 is over a very wide range dilution authentic, the linearity of serum dilutions is over a very wide range excellent (s.Tab.5).

**Table 5:** Recovery and linearity of the Sample Dilution (characteristic results of two different sera)

Dilution	Sample 1 (native 5.5 ng/ml)		Sample 2 (native 2.25 ng/ml)	
	plus 5 ng/ml	Recovery (%)	plus 12.25 ng/ml	Recovery (%)
1:50	9.71	92.5	14.99	103.4
1:100	10.60	101.0	13.64	94.1
1:200	10.44	99.4	14.10	97.2
1:400	10.32	98.3	14.33	98.8

Different human sera were spiked with recombinant human Resistin in varying concentrations (e.g. in Table 6). The recovery of Resistin yielded on average 98 % of the theoretically expected amount.

**Table 6:** Samples were enriched with 0.3 ng/ml Resistin and measured in comparison to non enriched sample. Relative recovery of added Resistin is shown.

Matrix	Dilution	% Recovery
Cerebrospinal fluid	1:2	129
Cerebrospinal fluid	1:10	93
Cerebrospinal fluid	1:40	103
Amnion fluid	1:10	85
Amnion fluid	1:40	91
Saliva	1:10	99
Saliva	1:21	86
Urine	1:10	79
Urine	1:21	85
Breast milk	1:2	97
Breast milk	1:10	58
Breast milk	1:21	63
Cell culture supernatant	1:2	100

## EVALUATION OF RESULTS

**Table 7:** The expected values for Resistin were determined with the Mediagnost ELISA E50 in healthy probands and analysed by Prof. Dr. J. Kratzsch, Institute for Laboratory Medicine, University of Leipzig.

Female				Resistin (ng/ml):		
Age (Years):	n:	AV Age:	AV BMI:	AV ± SD:	25.- 75. Percentile:	Min. – Max.:
18 - 30	96	23.0	23.1	7.2 ± 2.6	5.4 – 8.8	3.1 – 14.7
31 - 40	63	36.5	24.3	8.1 ± 2.3	6.4 – 9.6	3.6 – 13.1
41 - 50	67	44.9	24.8	7.3 ± 2.5	5.7 – 8.1	4.0 – 16.1
51 - 60	29	54.7	25.0	7.2 ± 2.6	5.4 – 8.5	4.0 – 15.5
61 - 65	9	62.7	25.2	6.6 ± 1.1	6.0 – 6.7	5.4 – 9.3
Male				Resistin (ng/ml):		
Age (Years):	n:	AV Age:	AV BMI:	AV ± SD:	25.- 75. Percentile:	Min. – Max.:
18 - 30	107	23.9	24.1	6.4 ± 1,8	5.0 – 7.6	2.5 – 13.1
31 - 40	59	35.9	25.0	6.7 ± 3,2	4.8 – 7.4	3.8 – 26.9
41 - 50	66	45.0	25.2	6.5 ± 2,8	4.5 – 7.4	2.4 – 16.7
51 - 60	36	54.8	26.4	6.1 ± 2,1	4.7 – 7.2	3.2 – 13.3
61 - 68	20	63.2	25.6	7.2 ± 1,8	6.0 – 8.2	4.5 – 11.2

n=Number of Proband, AV=Average Value, BMI=Body Mass Index (kg/m<sup>2</sup>), SD=Standard Deviation

**Table 8:** Grouping of the expected values

Sex	Number	Mean [ng/ml]	Standard deviation	2.5. Percentile	9.5. Percentile
Male	288	6.48	2.44	3.32	11.68
Female	264	7.41	2.47	3.68	13.60
Total	552	6.93	2.49	3.58	13.12

## LITERATUR / LITERATURE

1. Nakano, Y., et al., *Isolation and characterization of GBP28, a novel gelatin-binding protein purified from human plasma*. J Biochem (Tokyo), 1996. **120**(4): p. 803-12.
2. Pajvani, U.B., et al., *Structure-function studies of the adipocyte-secreted hormone Acrp30/Resistin. Implications for metabolic regulation and bioactivity*. J Biol Chem, 2003. **278**(11): p. 9073-85.
3. Tsao, T.S., et al., *Role of disulfide bonds in Acrp30/Resistin structure and signaling specificity. Different oligomers activate different signal transduction pathways*. J Biol Chem, 2003. **278**(50): p. 50810-17.
4. Shimada, K., T. Miyazaki, and H. Daida, *Resistin and atherosclerotic disease*. Clin Chim Acta, 2004. **344**(1-2): p. 1-12.
5. Higashiura, K., et al., *Correlations of Resistin level with insulin resistance and atherosclerosis in Japanese male populations*. Clin Endocrinol (Oxf), 2004. **61**(6): p. 753-9.
6. Spranger, J., et al., *Resistin is independently associated with insulin sensitivity in women with polycystic ovary syndrome*. Clin Endocrinol (Oxf), 2004. **61**(6): p. 738-46.
7. Zoico, E., et al., *Adipocytokines, fat distribution, and insulin resistance in elderly men and women*. J Gerontol A Biol Sci Med Sci, 2004. **59**(9): p. M935-9.
8. Ye, J.M., et al., *PPARalpha /gamma ragaglitazar eliminates fatty liver and enhances insulin action in fat-fed rats in the absence of hepatomegaly*. Am J Physiol Endocrinol Metab, 2003. **284**(3): p. E531-40.
9. Yamauchi, T., et al., *Resistin stimulates glucose utilization and fatty-acid oxidation by activating AMP-activated protein kinase*. Nat Med, 2002. **8**(11): p. 1288-95.
10. Blüher (Blueher), M., et al., *Total and high-molecular weight Resistin in relation to metabolic variables at baseline and in response to an exercise treatment program: comparative evaluation of three assays*. Diabetes Care, 2007. **30**(2): p. 280-5.
11. Delaigle, A.M., et al., *Induction of Resistin in skeletal muscle by inflammatory cytokines: in vivo and in vitro studies*. Endocrinology, 2004. **145**(12): p. 5589-97.
12. Winzer, C., et al., *Plasma Resistin, insulin sensitivity, and subclinical inflammation in women with prior gestational diabetes mellitus*. Diabetes Care, 2004. **27**(7): p. 1721-7.
13. Xydakis, A.M., et al., *Resistin, inflammation, and the expression of the metabolic syndrome in obese individuals: the impact of rapid weight loss through caloric restriction*. J Clin Endocrinol Metab, 2004. **89**(6): p. 2697-703.
14. Motoshima, H., et al., *Resistin suppresses proliferation and superoxide generation and enhances eNOS activity in endothelial cells treated with oxidized LDL*. Biochem Biophys Res Commun, 2004. **315**(2): p. 264-71.
15. Wolf, A.M., et al., *Resistin induces the anti-inflammatory cytokines IL-10 and IL-1RA in human leukocytes*. Biochem Biophys Res Commun, 2004. **323**(2): p. 630-5.
16. Okamoto, Y., et al., *Resistin reduces atherosclerosis in apolipoprotein E-deficient mice*. Circulation, 2002. **106**(22): p. 2767-70.
17. Schlegel, A., *Resistin and risk of coronary heart disease*. Jama, 2004. **292**(1): p. 40; author reply 40.
18. Choi, K.M., et al., *Inflammation, Insulin Resistance, and Glucose Intolerance in Acute Myocardial Infarction Patients without a Previous Diagnosis of Diabetes Mellitus*. J Clin Endocrinol Metab, 2004.
19. Nakamura, Y., et al., *Implications of plasma concentrations of Resistin in patients with coronary artery disease*. Heart, 2004. **90**(5): p. 528-33.
20. Pischon, T., et al., *Plasma Resistin levels and risk of myocardial infarction in men*. Jama, 2004. **291**(14): p. 1730-7.
21. Shibata, R., et al., *Resistin stimulates angiogenesis in response to tissue ischemia through stimulation of amp-activated protein kinase signaling*. J Biol Chem, 2004. **279**(27): p. 28670-4.
22. Fernandez-Real, J.M., et al., *Resistin is associated with vascular function independent of insulin sensitivity*. Diabetes Care, 2004. **27**(3): p. 739-45.

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## SUMMARY – MEDIAGNOST RESISTIN ELISA E50

Reagent preparation:	Reconstitution:	Dilution:
Standards A – E	in 750 µl Sample Buffer PP	
Control Serum KS	in 100 µl Dilution Buffer VP	1:21 with Sample Buffer PP
Antibody Conjugate AK		1:100 with Dilution Buffer VP
Enzyme Conjugate EK		1:100 with Dilution Buffer VP
Washing Buffer WP		1:20 with Aqua.dest. (e.g., add the complete contents of the flask (50 ml) into a graduated flask and fill with A.dest. to 1000 ml).
<b>Sample dilution: 1:21 (e.g. 15 µl Serum with 300 µl Sample Buffer PP)</b>		

### Assay Procedure for Double Determination

Pipette	Reagents	Position
100 µl	Sample Puffer PP (blank value)	A1/2
100 µl	Standard A (0.02 ng/ml)	B1/2
100 µl	Standard B (0.1 ng/ml)	C1/2
100 µl	Standard C (0.3 ng/ml)	D1/2
100 µl	Standard D (0.6 ng/ml)	E1/2
100 µl	Standard E (1.0 ng/ml)	F1/2
100 µl	Control Serum KS	G1/2
100 µl	Sample dilution	following wells
Cover the wells with the sealing tape.		
<b>Incubation: 2 h at RT, 350 rpm</b>		
5x 300 µl	Aspirate the contents of the wells and wash 5x with 300 µl Wash Buffer WP	each well
100 µl	1:100 diluted Antibody Conjugate AK	each well
<b>Incubation: 1 h at RT, 350 rpm</b>		
5x 300 µl	Aspirate the contents of the wells and wash 5x with 300 µl Wash Buffer WP	each well
100 µl	1:100 diluted Enzyme Conjugate EK	each well
<b>Incubation: 30 min at RT, 350 rpm</b>		
5x 300 µl	Aspirate the contents of the wells and wash 5x with 300 µl Wash Buffer WP	each well
100 µl	Substrate Solution S	each well
<b>Incubation: 30 min in the dark at RT</b>		
100 µl	Stop Solution SL	each well
Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength		

REF E50

## International Test description

CAL A-E	A -E	Rec in 750 µl	BUF PP	
Control	KS	Rec in 100 µl	BUF VP	1:21 DILU BUF PP
AB	AK			1:100 DILU BUF VP
CONJ	EK			1:100 DILU BUF VP
WASHBUF 20x	WP			1:20 DILU A. dest.
SPE				1:21 DILU BUF PP
°C 20-28 °C				

100 µl	BUF PP		A1/2
100 µl	CAL A (0.02 ng/ml)		B1/2
100 µl	CAL B (0.1 ng/ml)		C1/2
100 µl	CAL C (0.3 ng/ml)		D1/2
100 µl	CAL D (0.6 ng/ml)		E1/2
100 µl	CAL E (1 ng/ml)		F1/2
100 µl	CONTROL KS 1:21 DILU BUF PP		G1/2
100 µl	SPE 1:21 DILU BUF PP		
TAPE			

2 h 20-25 350 rpm

5x 300 µl	5x WASHBUF WP
100 µl	Ab AK
TAPE	

1 h 20-25 350 rpm

5x 300 µl	5x WASHBUF WP
100 µl	CONJ EK
TAPE	

0.5 h 20-25 350 rpm

5x 300 µl	5x WASHBUF WP
100 µl	SUBST TMB S

30 min 20-25

100 µl	H <sub>2</sub> SO <sub>4</sub> SL
MEASURE	