OSTEOPROTEGERIN

(EN) ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF
OSTEOPROTEGERIN IN EDTA PLASMA, HEPARIN PLASMA, SERUM OR CELL
CULTURE SUPERNATANTS
CAT. NO. BI-20402 . 12 X 8 TESTS

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www.bmgrp.com

1) INTRODUCTION

Osteoprotegerin (OPG) or Osteoclast inhibitory factor (OCIF) is a dimeric glycoprotein of the TNF receptor family with a molecular weight of 60 kD and 120 kD respectively. OPG acts as a soluble secreted receptor for RANKL. OPG prevents the binding of RANKL to RANK and thus inhibits RANK activation (osteoclast differentiation and activation receptor, ODAR) on the osteoclast surface. OPG is an inhibitor of osteoclast development.

Indications:

- Osteoporosis (1, 2)
- Diseases with locally increased resorption activity (3- 5)
- Arthritis (6)
- Oncology (7-10)
- Therapy monitoring (11, 12)
- Cardiovascular Disease (13-16)

2) CONTENT OF THE KIT

CONT	KIT COMPONENTS	QUANTITY
PLATE	Mouse monoclonal anti OPG antibody, pre-coated microtiter strips in a strip holder	12 x 8 tests
WASHBUF	Wash buffer concentrate 20x, natural cap	1 x 50 ml
AB	Goat polyclonal anti OPG antibody – biotin labelled, natural cap, yellow dye, ready to use	1 x 7 ml
STD	Standards 1-6, (0; 0.37; 1.1; 3.3; 10; 30 pmol/l), white caps, ready to use Only for use with plasma or serum samples	6 x 500 µl
CTRL	Control, yellow cap, ready to use (exact concentration on the label) Only for use with plasma or serum samples	1 x 500 µl
ASYBUF	Assay Buffer, red cap, ready to use	1 x 25 ml
CONJ	Conjugate, (streptavidin-HRPO), amber bottle, amber cap, ready to use	1 x 22 ml
SUB	Substrate (TMB solution), amber bottle, blue cap, ready to use	1 x 22 ml
STOP	Stop solution, white cap, ready to use	1 x 7 ml
OPG-STOCK	Osteoprotegerin stock solution, red cap, (500 pmol/l) Only for use with cell culture samples, (see chapter 5) reagents and sample preparation	1 x 500 µl

3) ADDITIONAL MATERIAL IN THE KIT

- 2 self-adhesive plastic films
- Quality control protocol
- Protocol sheet
- Instruction for use

4) MATERIAL AND EQUIPMENT REQUIRED BUT NOT SUPPLIED

- 1.5 ml reaction vials
- Precision pipettes calibrated to deliver 50 µl, 100 µl, 200 µl, 250 µl, 300 µl, 1000 µl, and disposable tips
- Distilled or deionised water
- Plate washer is recommended for washing, alternative multichannel pipette or manifold dispenser
- Refrigerator with 4°C (2-8°C)
- ELISA reader capable of measuring absorbance at 450 nm (with correction wavelength at 630 nm)
- Graph paper or software for calculation of results

5) REAGENTS AND SAMPLE PREPARATION

All reagents of the kit are stable at 4°C (2-8°C) until expiry date stated on the label of each reagent.

Sample preparation:

Collect venous blood samples by using standardized blood collection tubes for serum or plasma. We recommend performing plasma or serum separation by centrifugation as soon as possible (e.g. 20 min at 2000 x g, preferably at 4°C (2-8°C)). If this is not possible store the samples at 4°C (2-8°C) prior to centrifugation (up to one day). The acquired plasma or serum samples should be measured as soon as possible. For longer storage aliquot samples and store at -20°C or at -70°C. All samples should undergo only 4 freeze-thaw cycles. Lipemic or haemolysed samples may give erroneous results. Samples should be mixed well before assaying.

For further information on sample stability please visit our website www.bmgrp.com/products/assay_characteristics/sample_stability or contact our customer service by e-mail export@bmgrp.at or by phone +43/ 1/ 29107-45.

Reconstitute as follows:

WASHBUF (Wash buffer): Salt precipitate in the concentrated wash buffer is normal. Dissolve any precipitate by mixing gently at room temperature then dilute the concentrate 1:20 with distilled/DI water (e.g. 50 ml WASHBUF + 950 ml distilled water). prior to using in the assay. Diluted wash buffer is stable at 4°C (2-8°C) for one month. Only use diluted WASHBUF (Wash buffer) for optimum assay performance.

Measurement of Cell culture samples: OPG STOCK (OPG stock solution)

Only for measurement of OPG in Cell Culture Samples. Do not use serum standards 1-6 (white caps) and control (yellow cap)!

- Dilute the OPG-STOCK with cell culture medium 1:20 (1+19) for STD1 (50 µl OPG-STOCK + 950 µl medium)
- Dilute the standard for the standard curve 1:2 (1+1) as it is shown in the scheme:

STD1 (25 pmol/l)

250 µl STD1 + 250 µl cell culture medium = STD2 (12.5 pmol/l)

250 µl STD2 + 250 µl cell culture medium = STD3 (6.25 pmol/l)


250 µl STD3 + 250 µl cell culture medium = STD4 (3.13 pmol/l)

Use cell culture medium as STD0 (Standard 0)

Diluted OPG STOCK is stable at -20°C for one month. Avoid freeze-thaw cycles.

The range of standards depends from the quantity of produced OPG in the cell line. If necessary dilute the supernatant of the cell culture with cell culture medium.

6) PRINCIPLE OF THE ASSAY

7) ASSAY PROTOCOL


All reagents and samples must be at room temperature (18-26°C) before use in the assay
Mark position for BLANK/STD/SAMPLE/CTRL (Blank/Standard/Sample/Control) on the protocol sheet
Take microtiter strips out of the alu bag, take a minimum of one well as Blank. Store unused strips with desiccant at 4°C (2-8°C) in the alu bag. Strips are stable until expiry date stated on the label

- 1) Pipette 100 µl ASYBUF (Assay Buffer, red cap) into each well. Pipette additional 100 µl into well marked as blank
- 2) Add 50 µl STD/SAMPLE/CTRL (Standard/Sample/Control) in duplicate into respective well, except blank
- 3) Add 50 µl AB (biotinylated anti OPG antibody, natural cap) into each well, except blank, swirl gently
- 4) **Cover tightly and incubate over night (18-24h) at 4°C (2-8°C)**
- 5) Aspirate and wash wells 5x with 300 µl diluted WASHBUF (Wash buffer). After final wash, remove remaining WASHBUF by strongly tapping plate against paper towel
- 6) Add 200 µl CONJ (Conjugate, amber cap) into each well
- 7) **Cover tightly and incubate for 1 hour at room temperature (18-26°C)**
- 8) Aspirate and wash wells 5x with 300 µl diluted WASHBUF (Wash buffer). After final wash, remove remaining WASHBUF by strongly tapping plate against paper towel
- 9) Add 200 µl SUB (Substrate, blue cap) into each well
- 10) **Incubate for 20 min at room temperature (18-26°C) in the dark**
- 11) Add 50 µl STOP (Stop solution, white cap) into each well
- 12) Measure absorbance immediately at 450 nm with reference 630 nm, if available

8) CALCULATION OF RESULTS

Read the optical density (OD) of all wells on a plate reader using 450 nm wavelength (correction wavelength 630 nm). Subtract the blank OD from the values of STD, CTRL and sample. Construct the standard curve from the OD values of the STD. For cell culture samples use the standards prepared from the OPG-STOCK, see chapter 5) reagents and sample preparation. Use commercially available software or graph paper. Obtain sample concentration from this standard curve. The assay was evaluated with 4PL algorithm. Different curve fitting methods need to be evaluated by the user.

Example typical STD-curve:

The quality control (QC) protocol supplied with the kit shows the results of the final release QC for each kit at production date. Data for OD obtained by customers may differ due to various influences and/or due to the normal decrease of signal intensity during shelf life. However, this does not affect validity of results as long as an OD of 1.50 or more is obtained for the STD with the highest concentration and the value of the CTRL is in range (target range see label).

9) ASSAY CHARACTERISTICS

Normal range:	Median = 1.8 pmol/l (n = 1134; age 19-96) It is recommended to establish the normal range for each laboratory.
Standard range:	0 to 30 pmol/l
Conversion factor pg/ml to pmol/l:	1 pg/ml = 0.05 pmol/l (MW: 19.9 kD)
Sample volume:	50 µl human EDTA plasma, Heparin plasma, serum, cell culture supernatant
Detection Limit:	(0 pmol/l + 3 SD): 0.14 pmol/l
Incubation time:	18-24 h / 1 h / 20 min

For further information on assay characteristics please visit our website www.bmgrp.com/products/assay_characteristics or contact our customer service by e-mail export@bmgrp.at or by phone +43/ 1/ 29107-45.

10) PRECISION

Intra-Assay: 2 samples of known concentrations were tested 16 times to assess intra-assay precision

Inter-Assay: 2 samples of known concentrations were tested in 10 assays to assess inter-assay precision

Intra-Assay (n=16)	Sample 1	Sample 2
Mean (pmol/l)	4.59	10.76
SD	0.46	0.41
CV%	10%	4%

Inter-Assay (n=10)	Sample 1	Sample 2
Mean (pmol/l)	5.53	10.1
SD	0.38	0.76
CV%	7%	8%

11) TECHNICAL HINTS

- Do not mix or substitute reagents with those from other lots or sources.
- Do not mix stoppers and caps from different reagents or use reagents between lots.
- Do not use reagents beyond expiration date.
- Protect reagents from direct sunlight.
- Substrate solution should remain colourless until added to the plate.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.

12) PRECAUTIONS

All test components of human source were tested with 3rd generation tests against HIV-Ab and HBsAg, and were found negative. Nevertheless, they should be handled and disposed as if they were infectious.

All liquid reagents contain 0.01% Proclin 300 as preservative. Avoid contact with skin and mucous membrane.

Proclin 300 is not toxic in concentrations used in this kit. It may cause allergic skin reactions – avoid contact with skin or eyes.

- Do not pipette by mouth.
- Do not eat, drink, smoke or apply cosmetics where reagents are used.
- Wear gloves, glasses and lab jacket while performing this assay.
- Sulfuric acid is irritating to eyes and skin. Avoid contact with skin and mucous. Irritations are possible - Flush with water if contact occurs!!

13) LITERATURE

- 1) Samelson EJ et al.: Increased Plasma Osteoprotegerin Concentrations are Associated with Indices of Bone Strength of the Hip. *J Clin Endocrinol Metab* (2008), 93: 1789-1795.
- 2) Marini H et al.: OPG and sRANKL Serum Concentrations in Osteopenic, Postmenopausal Women after 2-year Genistein Administration. *JBMR* (2008), 23, 5: 715-720.
- 3) Kastritis E et al.: Increased Osteoprotegerin Levels in Response to Increased Bone Resorption in Amyloidosis: Comparisons with Bone Remodeling of Multiple Myeloma. *Blood* (2007), 110: 3528.
- 4) Madarász E et al.: Osteoprotegerin Levels in Women With Prior Gestational Diabetes Mellitus. *Diabetes Care* (2009), 32: e5.
- 5) Kearns AE et al.: Receptor Activator of Nuclear Factor κ B Ligand and Osteoprotegerin Regulation of Bone Remodeling in Health and Disease. *Endocr Rev* (2008), 29: 155 - 192.
- 6) Anastasilakis AD et al.: Evaluation of Bone Mineral Density, Bone Metabolism, Osteoprotegerin and Receptor Activator of the NF κ B Ligand Serum Levels During Treatment with Infliximab in Patients with Rheumatoid Arthritis. *Eur J Endocrinol* (2008), 158: 411 – 415
- 7) Corey E et al.: Osteoprotegerin in Prostate Cancer Bone Metastasis. *Cancer Res* (2005), 65: 1710-1718.
- 8) Zojer N et al.: Bisphosphonate Treatment does not Affect Serum Levels of Osteoprotegerin and RANKL in Hypercalcemic Cancer Patients. *Anticancer Res* (2005), 25: 3607-3612.
- 9) Terpos E et al.: The Clinical Significance of Serum Markers of Bone Turnover in NSCLC Patients: Surveillance, Management and Prognostic Implications.: *Anticancer Res* (2009), 29: 1651 - 1657.
- 10) Rogers A et al. (Eastele): Circulating Osteoprotegerin and Receptor Activator for Nuclear Factor κ B Ligand: Clinical Utility in Metabolic Bone Disease Assessment. *J Clin Endocrinol Metab* (2005), 90: 6323-6331.
- 11) Felicia Faienza M et al.: Osteoclastogenesis in Children with 21-Hydroxylase Deficiency on Long-Term Glucocorticoid Therapy: The Role of Receptor Activator of Nuclear Factor- κ B Ligand/Osteoprotegerin Imbalance. *J Clin Endocrinol Metab* (2009), 94: 2269-2276.
- 12) Woo JH et al.: Changes of Clinical Response and Bone Biochemical Markers in Patients with Ankylosing Spondylitis taking Etanercept. *J Rheumatol* (2007), 34: 1753 - 1759.
- 13) Semb A et al.: Osteoprotegerin and Soluble Receptor Activator of Nuclear Factor- κ B Ligand and Risk for Coronary Events. A Nested Case-Control Approach in the Prospective EPIC-Norfolk Population Study 1993-2003. *Arterioscler Thromb Vasc Biol* (2009), 29: 975-980.
- 14) Schroff RC et al.: The Circulating Calcification Inhibitors, Fetuin-A and Osteoprotegerin, but not Matrix Gla protein, are Associated with Vascular Stiffness and Calcification in Children on Dialysis. *Nephrol Dial Transplant* (2008), 23: 3263-3271
- 15) Morena M et al.: Cut-off value of Plasma Osteoprotegerin Level May Predict the Presence of Coronary Artery Calcifications in Chronic Kidney Disease Patients. *Nephrol. Dial. Transplant.* (2009), 10.1093/ndt/gfp301.
- 16) Omland T et al.: Circulating Osteoprotegerin Levels and Long-Term Prognosis in Patients With Acute Coronary Syndromes. *J Am Coll Cardiol* (2008), 51:627-633.

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ASSAY PROTOCOL AND CHECKLIST

PREPARATION OF REAGENTS:

- Bring all reagents to room temperature (18-26°C).
- Prepare reagents and samples as instructed.
- Bring unused and prepared components to the storage temperature mentioned in the package insert.
- Take microtiter strips out of the alu bag and mark positions on the protocol sheet.

TEST PROCEDURE:

- Step 1) Pipette 100 µl ASYBUF (Assay buffer, red cap) into each well. Pipette additional 100 µl into well marked as blank.
- Step 2) Add 50 µl STD/SAMPLE/CTRL (standard/sample/control) into all wells except blank.
- Step 3) Add 50 µl AB (biotinylated anti OPG ab, natural cap) into each well, except blank, swirl gently.
- Step 4) Cover tightly and incubate over night (18-24 hours) at 4°C.**
- Step 5) Aspirate and wash wells with 300 µl WASHBUF (Wash buffer) five times. Remove remaining buffer by hitting plate against paper towel.
- Step 6) Add 200 µl CONJ (Conjugate, amber cap) into each well.
- Step 7) Cover tightly and incubate for 1 hour at room temperature (18-26°C).**
- Step 8) Aspirate and wash wells with 300 µl WASHBUF (Wash buffer) five times. Remove remaining buffer by hitting plate against paper towel.
- Step 9) Add 200 µl SUB (Substrate, blue cap) into each well.
- Step 10) Incubate for 20 minutes at room temperature (18-26°C), in the dark.**
- Step 11) Add 50 µl STOP (Stop solution, white cap) into each well.
- Step 12) Read Optical Density at 450 nm with reference 630 nm, if available.