Ghrelin – RIA

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PACKAGE INSERT ENGLISH

TECHNICAL FEATURES+APPLICATIONS

- ◆ analytical sensitivity 0.04 ng/ml
- ◆ Intra- and Inter Assay Variance < 10%
- ◆ Recovery of recombinant Ghrelin 97%
- ◆ Control Serum included
- ◆ For measurement in human serum

INTENDED USE

This radioimmunoassay kit is suited for measuring human Ghrelin in serum and EDTA-Plasma.

INTRODUCTION

Ghrelin is a 3.5 kDa protein of 28 amino acids and is serine octanoylated. Bioactivity of this small peptide hormone depends on octanoylation (1). It is mainly synthesized by stomach but also in duodenal and heart cells (2) and therefore indicates the relevance of the stomach as endocrine organ (3).

Ghrelin is able to cross the blood-brain barrier and is a natural ligand of growth hormone secretagogue receptor in pituitary and hypothalamus (4). Ghrelin exert influence on several neurological processes, for instance memory retention can be modulated by ghrelin (5) and Ghrelin secretion is influenced by sleep (6, 7). Not only growth hormone but several other hormones are influenced by Ghrelin e.g. ACTH, cortisol, prolactin (8, 9). It is also present in pancreatic islets and regulates insulin secretion (10-12). In women with polycystic ovary syndrome Ghrelin levels are decreased and highly correlated to insulin sensitivity (24), so there are several regulatory circles influenced by Ghrelin. Additionally to endocrine action Ghrelin exerts influence on immunological processes. In Human umbilical vein endothelial cells Ghrelin inhibits basal and TNF-alpha-induced cytokine release and mononuclear cell binding and in vivo endotoxin-induced proinflammatory cytokine production in rats was also inhibited by intravenous administrated Ghrelin (13).

Sites of Ghrelin synthesis as well as receptor location indicate a role for the hormone and the gut in appetite regulation (1) (23). So many Ghrelin receptors are present in the hypothalamic arcuate nucleus, a brain area important in food intake control. Several investigations demonstrate a circadian rhythm of Ghrelin secretion (14), controlled by ingestion. Shortly before food intake Ghrelin plasma concentration increases and decreases after finishing. In eating disorders Ghrelin levels reflect illness, so obesity suppresses Ghrelin concentration in blood (14) and in anorexia nervosa an increase of Ghrelin serum concentration can be detected (15-17). Ghrelin might act as counterpart to Leptin in the regulation of food intake and fat utilization, so in patients with primary biliary cirrhosis parallel to increasing Leptin levels Ghrelin serum concentration decreases (16). It also influences the adipogenesis negatively (18). A significant decrease of Ghrelin concentration is detected in elderly people (19). This could explain the anorexia of elderly people and offers a new target in anti-aging research.

Further investigation of Ghrelin serum concentration in chronic liver disease resulted in a correlation of Ghrelin concentration with several biochemical and clinical parameters: e.g. encephalopathy, anaemia, hypoglycaemia, renal dysfunction (20).

In Patients with the Prader-Willi-Syndrome show a nearly 5fold increase of serum Ghrelin, so Ghrelin may be involved in the pathogenesis of hyperphagia in PWS (21, 22).

So Ghrelin seems to be involved in the regulation of many physiological processes and its influence on many of these processes has not been investigated in detail. With this test system we provide an easy but reliable tool for the quantification of human Ghrelin in serum.

PRINCIPLE

For the Radioimmunoassay for the determination of human Ghrelin a polyclonal rabbit-antibody of high specificity is used. Ghrelin is measured quantitatively.

Standards are prepared from recombinant Ghrelin, ¹²⁵I-Tracer from a C-terminated peptide (AS 15 – 28) with presynthesised Tyrosine.

Calibration of the Assay

The assay was calibrated against the internal test of Medical School Hannover, Prof. Brabant.

PERFORMANCE CHARACTERISTICS

Sensitivity

The analytical **sensitivity** of the assay yields **40 pg/ml (about 10 pmol/L)** measured as 2x SD of zero standards.

Specificity

This assay is specific for human Ghrelin. The antibody shows in addition cross-reactivity with: rabbit, cat, guinea pig, hamster, goat, sheep, rat, horse, donkey, pig, dog, rabbit, mouse and bovine. No cross-reactivity was found with other proteins such as insulin or GH.

Table 1: Intra-Assay-Variation

	Number of determinations	Mean value [pg/ml]	VC%
Sample 1	6	1028.13	2.6
Sample 2	6	1275.95	4.0
Sample 3	6	1325	5.3

Table 2: Inter-Assay-Variation

	Mean value (pg/ml)	Standard deviation	VC%
Sample 1	762	55.3	7.3
Sample 2	1085	52.5	4.8
Sample 3	790	65	8.2

Recovery

Serum spiking experiments with recombinant human Ghrelin yielded a recovery of 97% (± 2%).

SPECIMEN COLLECTION, PREPARATION AND STORAGE

Serum/ EDTA-Plasma samples are suitable (inappropriate are Heparin- and Citrate-Plasma). An external sample preparation prior to assay is not required.

For testing of single blood samples, the specimens may be taken in the morning or early afternoon. For specific questions, the influence of food intake should be taken in consideration.

Samples should be handled as recommended in general: as fast as possible and chilled as soon as possible. In case there will be a longer period between the sample withdrawal and determination store the undiluted samples frozen at –20°C or below in tightly closable plastic

tubes. Avoid on principal repeated freeze-thaw cycles of serum/plasma (if required, please subaliquote) although Ghrelin levels were found to be unaffected by few cycles (5x) in our experiments.

Because of the wide effective range of this RIA kit a preparative sample dilution is generally not necessary. For most of the determinations (serum or plasma samples, and no extreme values expected) **the use of undiluted samples 100 µl per tube**, should be appropriate. In case of extremely high Ghrelin levels the sample should and can be diluted in Assay Buffer AB, e.g. 1:10. At a concentration of 1000pg/ml the precision of this assay is maximal.

REAGENTS PROVIDED

1)	CAL 1-6	Standards STD 1-6 , lyophilized, 750 µl , contain human recombinant Ghrelin. The calibration curve covers a range of 2-6.4 ng/ml Ghrelin. Please use 100 µl standard solution per tube .
2)	NSB	unspecific binding, NSB , lyophilized, 1ml . rabbit IgG
3)	BUF AB	Assay Buffer AB , ready-to-use , 30 ml , for reconstitution of 1.AB, TR, NSB, STD and CS
4)	BUF PR	Precipitation-Reagent PR , after addition of 2.Ab ready-to use, dilution of 2.Ab 1:56, 55ml
5)	1.Ab	1. antibody, 1.Ab , lyophilized, 10,5 ml , Reconstitute in 10.5 ml AB.
6)	2.Ab	2. antibody, 2.Ab , lyophilized, 1 ml , anti-rabbit-IgG. Reconstitute in 1 ml AB. 2.Ab Reconstitute with 1 ml reagent AB . Transfer dissolved material to reagent PR immediately before use (Ratio 1:56). Please mix only the required quantity of PR+2.Ab (500µl/vial). The rest can be frozen. The assay is unaffected by the possible occurrence of turbidity after adding 2.Ab to reagent PR .
7)	Tracer	Tracer TR , lyophilized, 10.5 ml , < 2.3 µCi or < 85 kBq. reconstitute in 10.5 ml AB, stained red
8)	Control	Control CS , lyophilized, 750 µl , reconstitute in 750 µl AB. Concentration is given by label +/- 2SD.

MATERIALS REQUIRED BUT NOT PROVIDED

Precision pipettes (100 and 200µl) Micropipettes and multichannel pipettes with disposable plastic tips
 Disposable polystyrene or polypropylene tubes. Conical tubes are highly recommended because of the small volume of the immunoprecipitates

Vortex-mixer

Centrifuge

Device to aspirate the fluid from the tubes (recommended because of the potential danger of radioactivity and infection by human samples)

Ice-Cold deionized water

Gamma Counter

REAGENT PREPARATION

In conducting the assay, follow strictly the test protocol. Room temperature incubation means: Incubation at 20 - 25°C.

Reagents with different lot numbers should not be mixed. All reagents are stable unopened until the expiry date, if stored in the dark at 2° - 8°C (see label).

Control Serum CS Reagent 1.Ab, 2.Ab, TR, NSB and STD have to be reconstituted in **Assay Buffer AB**. It is recommended to keep the reconstituted reagents at room temperature for 30 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer. Transfer reconstituted 2.Ab into the buffer PR. The assay is unaffected by the possible occurrence of turbidity after add in 2.Ab to reagent PR.

The shelf life of the components after opening is not affected, if used appropriately. Reconstituted Components should be stored at -20°C (or below). Repeated freeze-thaw cycles have to be avoided.

Before use, all kit components should be brought to room temperature, if nothing different is indicated. **Precipitates, possible in buffers, should be dissolved before use through mixing and warming.**

WARNINGS AND PRECAUTIONS

For in-vitro diagnostic use only. For professional use only.

Possession and use of the kit is subject to the regulations of the national nuclear regulatory authorities.

Reagents with different lot numbers should not be mixed.

Reagents contain Sodium-Azide as preservative, however, highly diluted (0.02%). Sodium-Azide is very toxic, R-Phrases: 28, 32, 50/53 and S-Phrases 28, 45, 60, 61 must be considered.

Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.

Before use, all kit components should be brought to **room temperature at 20 - 25°C**, if not indicated differently. Precipitates in buffers should be dissolved before use by thorough mixing and warming. **Temperature WILL affect** the assay. However, values for the patient samples will not be affected.

Caution: This kit contains material of human and/or animal origin. Source human serum for the Control Serum provided in this kit was tested by FDA recommended methods and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV) antibodies. No known test methods can offer total assurance of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

Radioactivity - Before ordering or using radioactive materials, it is necessary to take the appropriate actions to ensure compliance with national regulations governing their use. Local rules in each establishment, which define actions and behaviour in the radioactivity working areas, should also be adhered to. The advice given here does not replace any local rules, instructions or training in the establishment, or advice from the radiation protection advisers. It is important to follow the code of good laboratory practice in addition to the specific precautions relating to the radionuclide I-125 used.

Iodine-125 has a radioactive half-life T_{1/2} of 60 days and emits 35.5 keV gamma radiation, 27 – 32 keV x-rays and no beta radiation. Shielding is effectively done by lead, first half value layer is 0.02 mm lead, reduction to 10 % is made by 0.2 mm.

To reduce the radiation dose time spent handling radioactivity should be minimized (plan ahead), and distance from source of radiation should be maximized (doubling the distance from the source quarters the radiation dose).

Formation of aerosols, e.g. by improper opening and mixing of vials or pipetting of solutions which may cause minute droplets of radioactivity become airborne, is a hazard and should be avoided.

Solutions containing iodine should not be made acidic, because this might lead to the formation of volatile elemental iodine.

As some iodo-compounds can penetrate rubber gloves, it is advisable to wear two pairs, or polyethylene gloves over rubber.

For cleaning of contaminated areas or equipment, the Iodine-125 should be rendered chemically stable by using alkaline sodium thiosulphate solution together with paper or cellulose tissue.

General first aid procedures:

Skin contact: Wash affected area thoroughly with water. Discard contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: If swallowed, wash out mouth thoroughly with water. Immediately see a physician.

Do not eat, drink or smoke in these areas.

Never pipette the materials with the mouth.

Spilled material must be wiped off immediately and should become disinfected. Clean contaminated areas and equipment with a suitable detergent.

The handling of radioactive and potentially infectious material must comply with the following guidelines:

The material should be stored and used in a special designated area.

Do not eat, drink or smoke in these areas.

Never pipette the materials with the mouth.

Avoid direct contact with these materials by wearing laboratory coats and disposable gloves.

Spilled material must be wiped off immediately. Clean contaminated areas and equipment with a suitable detergent.

Unused radioactive material and radioactive waste should be disposed according to the recommendations of the national regulatory authorities.

ASSAY PROCEDURE

NOTES: All determinations (Standards, Control Serum and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

When performing the assay, the Standards, Control Serum and the samples should be pipette as fast as possible.

Flow Chart of Assay Protocol

Tube Nr. :	Contents	AB	STD, CS, Samples	NSB	1.Ab	TR	2.Ab in PR
1,2	Total Counts	-	-	-	-	100	-
3,4	NSB	100	-	100	-	100	500
5,6	B ₀ (zero standard)	100	-	-	100	100	500
7-18	Standards 1-6	-	100 STD1-6	-	100	100	500
19,20	Control	-	100 CS	-	100	100	500
21,22	Sample 1	-	100	-	100	100	500
23,24	Sample 2	-	100	-	100	100	500
etc.							

All volumes are given in µl.

- 1) Labelling of the assay tubes should be done in the following order (duplicates):
 - 1, 2 total counts (**TC**),
 - 3, 4 non specific binding (**NSB**)
 - 5, 6 zero standard (**B₀**),
 - 7-18 **Standards** 1 to 6
 - 19, 20 Control **CS**
 - 21,22 etc. **samples**.
- 2) Add **100 µl** of reagent **AB** (Assay Buffer) to tubes 3, 4 and 5,6.
- 3) Add **100 µl** of **Standards 1-6** to tubes 7-18:
 - 7, 8 Standard 1 (200 pg/ml)
 - 9, 10 Standard 2 (400 pg /ml); etc. up to 18.
- 4) Add **100 µl** Control **CS** to tubes 19, 20.
- 5) Add **100 µl** of **sample** to tubes 21, 22, etc.
- 6) Add **100 µl NSB** to tubes 3 and 4.
- 7) Add **100 µl** of Reagent **1.Ab** (1st. Antibody), beginning with tube 5.
- 8) Mix tubes with a Vortex-Mixer and incubate **overnight** at 2-8°C (**at least 20 h, maximal 24 h**).
- 9) Add **100 µl** of Reagent **TR** (Tracer) to all tubes.
Seal tubes 1 and 2 (total counts) with a stopper and remove until step 16.
- 10) Mix the remaining tubes with a Vortex-Mixer and **incubate overnight (at least 16 h, maximal 20 h) at 2-8°C**.
- 11) Add from Reagent **2.Ab** to **PR** (1:56, mix); **the reagent mix must be cooled (2-8 °C)**. Add **500 µl** Reagent **PR** (incl. 2.Ab), beginning with tube 3.
- 12) Mix tubes with a Vortex-Mixer and incubate for precipitation for **1 h at 2-8°C**.
- 13) Add **1 ml ice-cold water**.
- 14) Centrifuge at 2-4°C at **3000 x g for 20 min**.
- 15) Aspirate the supernatant. In order not to destroy or aspirate the small precipitate a rest of approx. 2 mm supernatant should be left over the precipitate. (Tip: Add limit stop to the aspirate needle). Depending on laboratory equipments and common laboratory practice supernatant can also be decanted carefully.
- 16) Count the radioactivity of all tubes in Gamma-Counter for **1 to 3 min**.

Extended washing procedure for increased precision:

After the Incubation of 1 hour (step 12) centrifuge the tubes (s. step 14) and aspirate the supernatant (s. step 15). Add directly 1 ml ice-cold water. This should not be done too vigorously

in order to keep the precipitate intact. **Do not mix again!** Centrifuge the tubes once again at 3000g for 5 min, aspirate the supernatant and count the radioactivity of all tubes in the gamma-counter.

This extended procedure results in a somewhat higher precision bound up with higher work expenditure. The higher precision may be relevant only in special cases.

CALCULATION OF RESULTS

Establishing the Standard Curve

Standard	B ₀	1	2	3	4	5	6
ng/ml	0	0.2	0.4	0.8	1.6	3.2	6.4
pg/ml	0	200	400	800	1600	3200	6400
pmol/l	0	59	119	237	475	949	1899

1. Calculate the average counts of each pair of tubes.
2. Subtract the average of NSB (NSB, Tubes 3 and 4) from the mean counts of the Standards, controls and samples. This corresponds to the corrected B values.
3. The corrected value of the zero standard (tubes 5 and 6) equals B₀.
4. Calculate the percent bound (%B/B₀): $\%B/B_0 = B/B_0 \times 100 \%$
5. Plot %B/B₀ versus the standard concentrations on a semi-logarithmic or logit-log paper respectively or per computer analysis.
6. For quality control calculate the percentage of %NSB/TC:
 NSB / TC (average counts of tubes 3 and 4 / average counts of tubes 1 and 2) $\times 100 \%$. It should be: %NSB/TC < 5%
7. Quality control, Calculate the %B₀/TC:
 B_0 (see step 3.) / TC (total counts) $\times 100 \%$.
 It should be: %B₀/TC > 20%

Evaluation of sample concentrations:

Read the Ghrelin concentration value (abscissa) corresponding to the %B/B₀-Wert of the sample as in the example given below:

Example:

Average counts of the tubes 3 and 4 (NSB): 482 cpm

Average counts of the tubes 5 and 6 (Zero standard, B₀): 8927 cpm

Average counts of the tubes 19 and 20 (Control CS): 4794 cpm

$$\%B/B_0 = \frac{\text{Sample} - \text{NSB}}{B_0 - \text{NSB}} \times 100\% = \frac{4794 - 482}{8927 - 482} \times 100\% = 51.1\%$$

For a 51.1% value on the y-axis (ordinate) corresponds in this example abscissa value of 1055 pg/ml.

Multiplication of this value determined graphically or by aid of a computer program with the dilution factor gives the ghrelin concentration of the sample.

If it is preferred to express the results as pmol/l, the values given as pg/ml have to be divided by 3.371 to obtain pmol/l.

Example: 1055 pg/ml : 3.371 = 313 pmol/l.

Concentration of the control:

The control CS should fit with the labelled concentration range. The measured Ghrelin concentration are only valid, if the measured value of the control is in the labelled concentration range. Otherwise process analysis is required and depending on the results the measured values are accepted or not.

EXPECTED VALUES

Concentration of Ghrelin in human sera varied from minimal 300 pg/ml (90 pmol/L) to maximal 4000 pg/ml (1200 pmol/L) so far. Whereas in most of the samples the concentration ranged between 600 pg/ml (180 pmol/L) and 1400 pg/ml (420 pmol/L).

52 sera or plasma from each healthy male and female blood donors, age 20 to 65 years, were measured regarding their Ghrelin concentration with Mediagnost Ghrelin RIA R90. However, there is no information available regarding the nutritional status of the blood donors. The results in pg/ml are presented in the table 3.

Therefore valid reference values should be generated by an well designed study and the values presented here can only be taken as initial estimation.

Table 3: Ghrelin serum concentration of male and female blood donors.

	Male	Female
n	52	52
Mean	801.3	1141.9
Standard Deviation	179.1	570.9
Min	511.0	454.4
Max	1,555.2	2,565.7
Median	770.5	1,032.3

LITERATUR/LITERATURE

1. Kojima M, Hosoda H, Date Y, Nakazato M, Matsuo H, Kangawa K. Ghrelin is a growth-hormone-releasing acylated peptide from stomach. *Nature* 1999;402(6762):656-60.
2. Date Y, Kojima M, Hosoda H, et al. Ghrelin, a novel growth hormone-releasing acylated peptide, is synthesized in a distinct endocrine cell type in the gastrointestinal tracts of rats and humans. *Endocrinology* 2000;141(11):4255-61.
3. Inui A, Asakawa A, Bowers CY, et al. Ghrelin, appetite, and gastric motility: the emerging role of the stomach as an endocrine organ. *Faseb J* 2004;18(3):439-56.
4. Tschop M, Smiley DL, Heiman ML. Ghrelin induces adiposity in rodents. *Nature* 2000;407(6806):908-13.
5. Carlini VP, Varas MM, Cragnolini AB, Schiöth HB, Scimonelli TN, de Barioglio SR. Differential role of the hippocampus, amygdala, and dorsal raphe nucleus in regulating feeding, memory, and anxiety-like behavioral responses to ghrelin. *Biochem Biophys Res Commun* 2004;313(3):635-41.
6. Dzaja A, Dalal MA, Himmerich H, Uhr M, Pollmacher T, Schuld A. Sleep enhances nocturnal plasma ghrelin levels in healthy subjects. *Am J Physiol Endocrinol Metab* 2004.
7. Van Cauter E, Latta F, Nedeltcheva A, et al. Reciprocal interactions between the GH axis and sleep. *Growth Horm IGF Res* 2004;14 Suppl A:10-7.
8. Takaya K, Ariyasu H, Kanamoto N, et al. Ghrelin strongly stimulates growth hormone release in humans. *J Clin Endocrinol Metab* 2000;85(12):4908-11.
9. Broglio F, Arvat E, Benso A, et al. Ghrelin: endocrine and non-endocrine actions. *J Pediatr Endocrinol Metab* 2002;15 Suppl 5:1219-27.
10. Poykko SM, Kellokoski E, Horkko S, Kauma H, Kesaniemi YA, Ukkola O. Low plasma ghrelin is associated with insulin resistance, hypertension, and the prevalence of type 2 diabetes. *Diabetes* 2003;52(10):2546-53.
11. Anderwald C, Brabant G, Bernroider E, et al. Insulin-dependent modulation of plasma ghrelin and leptin concentrations is less pronounced in type 2 diabetic patients. *Diabetes* 2003;52(7):1792-8.
12. Soriano-Guillen L, Barrios V, Lechuga-Sancho A, Chowen JA, Argente J. Response of circulating ghrelin levels to insulin therapy in children with newly diagnosed type 1 diabetes mellitus. *Pediatr Res* 2004;55(5):830-5.
13. Li WG, Gavrilu D, Liu X, et al. Ghrelin inhibits proinflammatory responses and nuclear factor-kappaB activation in human endothelial cells. *Circulation* 2004;109(18):2221-6.
14. Cummings DE, Purnell JQ, Frayo RS, Schmidova K, Wisse BE, Weigle DS. A preprandial rise in plasma ghrelin levels suggests a role in meal initiation in humans. *Diabetes* 2001;50(8):1714-9.
15. Otto B, Cuntz U, Fruehauf E, et al. Weight gain decreases elevated plasma ghrelin concentrations of patients with anorexia nervosa. *Eur J Endocrinol* 2001;145(5):669-73.
16. Tolle V, Kadem M, Bluet-Pajot MT, et al. Balance in ghrelin and leptin plasma levels in anorexia nervosa patients and constitutionally thin women. *J Clin Endocrinol Metab* 2003;88(1):109-16.
17. Soriano-Guillen L, Barrios V, Campos-Barros A, Argente J. Ghrelin levels in obesity and anorexia nervosa: effect of weight reduction or recuperation. *J Pediatr* 2004;144(1):36-42.
18. Zhang W, Zhao L, Lin TR, et al. Inhibition of adipogenesis by ghrelin. *Mol Biol Cell* 2004;15(5):2484-91.
19. Rigamonti AE, Pincelli AI, Corra B, et al. Plasma ghrelin concentrations in elderly subjects: comparison with anorexic and obese patients. *J Endocrinol* 2002;175(1):R1-5.
20. Tacke F, Brabant G, Kruck E, et al. Ghrelin in chronic liver disease. *J Hepatol* 2003;38(4):447-54.
21. Cummings DE, Clement K, Purnell JQ, et al. Elevated plasma ghrelin levels in Prader Willi syndrome. *Nat Med* 2002;8(7):643-4.
22. Haqq AM, Farooqi IS, O'Rahilly S, et al. Serum ghrelin levels are inversely correlated with body mass index, age, and insulin concentrations in normal children and are markedly increased in Prader-Willi syndrome. *J Clin Endocrinol Metab* 2003;88(1):174-8.
23. Hanusch-Enserer U, Brabant G, Roden M. Ghrelin concentrations in morbidly obese patients after adjustable gastric banding. *N Engl J Med*. 2003 May 22; 348(21):2159-60.
24. Schofl C, Horn R, Schill T, Schlosser HW, Müller MJ, Brabant G. Circulating Ghrelin levels in patients with polycystic ovary syndrome. *J Clin Endocrinol Metab*. 2002 Oct; 87(10): 4607-10

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SUMMARY – Mediagnost Ghrelin R90

Reconstitution of the Reagents		
1st Antibody (1.Ab)	in Assay Buffer (AB)	10.5 ml
Tracer (TR)	in Assay Buffer (AB)	10.5 ml
NSB	in Assay Buffer (AB)	1 ml
Standards (1-6)	in Assay Buffer (AB)	750 µl
Control serum CS	In Assay Buffer (AB)	750 µl
2nd Antibody (2.Ab)	in Assay Buffer (AB) Mix solution with PR . Mix only the required quantity of PR+2ndAb	1 ml (1 ml 2.Ab+ 55 ml PR) or Ratio (1:56)

Serum or plasma **samples** can be used **undiluted**.

Double Determinations		Addition of Reagents [µl]			
Tube Nr. :	Contents	AB	STD (1-6), CS, Samples	NSB	1.AB (1 st Antibody)
1,2	Total Counts	-	-	-	-
3,4	NSB	100	-	100	-
5,6	B ₀ (zero standard)	100	-	-	100
7-18	Standards (1-6)	-	100 STD (1-6)	-	100
19,20	Control serum (CS)	-	100 CS	-	100
21,22	Sample 1	-	100	-	100
23,24	Sample 2	-	100	-	100
etc.					

Mix all tubes with a Vortex-mixer.

Incubation overnight (at least 20 h, maximal 24 h) **at 2-8°C**

Add **100 µl** Reagent **TR** (Tracer) to all tubes.

Seal the tubes 1 and 2 (total counts) with a stopper and **remove until the step 16**.

Mix the remaining tubes with a Vortex-Mixer.

Incubation overnight (at least 16 h, maximal 20 h) **at 2-8°C**

Add from the Reagent **2.Ab to PR (1:56, mix)**, the reagent mix must be cooled (**2-8°C**)
Add **500 µl** reagent **PR** (inc. 2.Ab.), beginning with the tube 3.

Mix tubes with a Vortex mixer.

Incubation 1 h at **2-8°C**

Add **1 ml ice-cold water**, beginning with the tube 3.

Centrifugation at **2-4°C**, **□3000 x g**, **20 min**

Aspirate the supernatant
(a rest of approx. 2 mm of supernatant should be left over the intact precipitate).

Count the radioactivity of all the tubes in a Gamma-Counter.

CAL 1-6	STD 1-6	Rec in 750 µl AB	
Control	CS	Rec in 750 µl AB	
NSB	NSB	Rec in 1 ml AB	
2.Ab	2.Ab	Rec in 1 ml AB	DILU PR 1:56 °C 2- 8
1.Ab	1. Ab	Rec in 10.5 ml AB	
Tracer	TR	Rec in 10.5 ml AB	
SPE			100 µl

°C 20-25 °C

Tubes			BUF AB	CAL	NSB	1.Ab 1.Ab
1/2	Tracer TR = TC	-	-	-	-	-
3/4	NSB	-	100 µl	-	100 µl	-
5/6	B₀	-	100 µl	-	-	100 µl
7/8	CAL STD 1 (0,2 ng/ml)	-	-	100 µl	-	100 µl
9/10	CAL STD 2 (0,4 ng/ml)	-	-	100 µl	-	100 µl
11/12	CAL STD 3 (0,8 ng/ml)	-	-	100 µl	-	100 µl
13/14	CAL STD 4 (1,6 ng/ml)	-	-	100 µl	-	100 µl
15/16	CAL STD 5 (3,2 ng/ml)	-	-	100 µl	-	100 µl
17/18	CAL STD 6 (6,4 ng/ml)	-	-	100 µl	-	100 µl
19/20	CONTROL KS	100 µl	-	-	-	100 µl
21/22	SPE	100 µl	-	-	-	100 µl
↔						

🕒 20-24h °C 2- 8

100µl **Tracer** TR ↔

🕒 16-20h °C 2- 8

500µl **2.Ab** **DILU** PR 1:56 °C 2- 8 **ALL** **NOT** **TC** ↔

🕒 1 h °C 2- 8

1 ml A. dest °C 2- 8 **ALL** **NOT** **TC**

🌀 ≥ 3000× g 20 min °C 2-8

ASP **ALL** **NOT** **TC****MEASURE**