25(OH)-Vitamin D *direct* ELISA Kit

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Hersteller

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
25(OH)-Vitamin D *direct* ELISA Kit

*For the determination of 25(OH)-Vitamin D
in human serum*

Valid from 28.03.2011

K 2109

EU patent #EP1097132

Australian patent # 763458

Content

1. INTENDED USE	19
2. SUMMARY AND EXPLANATION OF THE TEST	19
3. MATERIAL SUPPLIED	20
4. MATERIAL REQUIRED BUT NOT SUPPLIED	20
5. PREPARATION AND STORAGE OF REAGENTS	21
6. SPECIMEN COLLECTION AND PREPARATION	21
7. ASSAY PROCEDURE	22
<i>Principle of the test</i>	22
<i>Test procedure</i>	22
8. RESULTS	24
9. LIMITATIONS	24
10. QUALITY CONTROL	24
<i>Expected values</i>	24
11. PERFORMANCE CHARACTERISTICS	26
<i>Precision and reproducibility</i>	26
<i>Specificity</i>	26
<i>Day-to-Day Variation</i>	27
<i>Dilution Recovery</i>	27
<i>Spiking Recovery</i>	28
<i>Correlation data</i>	29
12. PRECAUTIONS	29
13. TECHNICAL HINTS	30
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	30
15. REFERENCES	31

1. INTENDED USE

The Immundiagnostik ELISA is intended for the quantitative determination of the 25-OH-Vitamin D in serum and fresh plasma. For *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Vitamin D is a steroid hormone involved in the intestinal absorption of calcium and the regulation of calcium homeostasis. There are two different forms of vitamin D, named D₃ and D₂, which are very similar in structure. The D₂ is a synthetic product, which is predominantly absorbed by fortified food.

Physiological vitamin D₃ levels result not only from dietary uptake but can also be produced from a cholesterol precursor, 7-dehydrocholesterol, in the skin during sun exposure. In the liver, the vitamin is hydroxylated to 25-hydroxyvitamin D (25(OH)-vitamin D), the major circulating metabolite of vitamin D. Although 1,25-(OH)₂ vitamin D portrays the biological active form of vitamin D, which is synthesized in the kidney, it is widely accepted that the measurement of circulating 25(OH)-vitamin D provides better information with respect to patients vitamin D status and allows its use in diagnose hypovitaminosis (1,2).

Preanalytical stability of 25(OH)-vitamin D₃ in human blood or serum at room temperature: solid as a rock (Wielders and Wijnberg, 2009).

The concentration of 25(OH)-vitamin D decreases with age and a deficiency is common among elderly persons.

Clinical applications of 25(OH)-vitamin D measurements are the diagnosis and therapy control of postmenopausal osteoporosis, rickets, osteomalacia, renal osteodystrophy, pregnancy, neonatal hypocalcemia and hyperparathyroidism. In addition, a prevalence of subclinical vitamin D deficiency has been discussed in different European countries.

Vitamin D intoxication mostly occurs during a large intake of pharmaceutical preparations of Vitamin D and may lead to hypercalcemia, hypercalcuria and nephrocalcinosis in susceptible infants.

3. MATERIAL SUPPLIED

Catalogue No.	Content	Kit Components	Quantity
K 2109MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 2109WP	WASHBUF	ELISA wash concentrate 20x	50 mL
K 2109RS	RECSOL	Reconstitution solution	2 x 20 mL
K 2109RR	RELREAG	Releasing reagent	2 x 1 vial
K 2109AK	AB	Anti 25(OH)-vitamin D antibody ready to use	18 mL
K 2109ST	STD	Standards, ready to use	6 vials 300 µL each
K 2109KO	CTRL	Controls, ready for use (see specification for range)	2 vials 300 µL each
K 2109K	CONJ	Conjugate, peroxidase labeled, ready to use	22 mL
K 2109TMB	SUB	TMB substrate (Tetramethylbenzidine)	2 x 15 mL
K 2109 SD	SAMDIL	Sample Dilution Buffer	60 mL
K 2109AC	STOP	ELISA stop solution, ready to use	1 x 15 mL
K 2109FOL	FOL	Foil to cover the microtiter plate	2 x 1

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Deep freezer -20 °C
- Precision pipettors calibrated to deliver 10–1000 µl
- Horizontal microtiter plate shaker
- Multi-channel dispenser or repeating dispenser
- Vortex-Mixer
- Water bath or heating block
- Standard laboratory glass or plastic vials, cups, etc. (one time products) made of polypropylene
- Microtiter plate reader 450 nm (reference wave length 620 or 690 nm)

20 ▪ Refrigerator with **defined 8–10 °C**

5. PREPARATION AND STORAGE OF REAGENTS

- The test kit is designed for 96 single determinations. A reduction of the sample or buffer volumes results in erroneous values.
- To run assay more than once, ensure that reagents are stored at conditions stated on the label. The kit can be used up to 2 times within the expiry date stated on the label.
- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:20** before use (50 ml WASHBUF + 950 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or at 37°C in a water bath before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. **Diluted buffer solution** can be stored in a closed flask at **2-8°C for two weeks**.
- **Pre-heat RECSOL (reconstitution solution) for at least 20 min at 37°C in a water bath before use.**
- Reconstitute **RELREAG** (releasing reagent) in **16 ml pre-heated RECSOL** (reconstitution solution), mix gently by careful swinging (do not vortex). The reconstituted RELREAG can be directly added into the prepared samples or kept at 37 °C in a water bath until use. After use, aliquot, freeze remaining releasing reagent and store at -20 °C. Frozen RELREAG can be defrosted and used only once. Pre-heat frozen releasing reagent to 37° C before use (e.g. incubate for **at least 20 minutes at 37° C in a waterbath**). Subsequently, it can be used right away.
- Bring **AB** (antibody) at room temperature at least one hour before use.
- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at **2 -8 °C**. **Note:** 1. Unused strips are covered and stored at 2-8° C. The covered strips should be used within 4 weeks. 2. The remaining RELREAG is stored at -20 °C.

6. SPECIMEN COLLECTION AND PREPARATION

1. Fresh collected blood should be centrifuged within one hour. Vitamin D is an inert substance. However, serum storage at 2-8°C is recommended when the analysis is performed within 24 h after collection. Otherwise, the serum samples must be stored at -20°C until analyzed. Avoid repeated freeze-thaw cycles.
2. Serum samples can be shipped at 4-8 °C (for example with Coolpacks) and remain stable for up to 3 days.
3. **Serum** is the preferred sample matrix; whole blood is not suitable.
4. Indicated incubation times and temperatures must be strictly observed.
5. Mix samples well before use.

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes of a competitive ELISA technique with a selected monoclonal antibody recognizing 25(OH)-vitamin D. For a reliable determination of 25(OH)-vitamin D, it is necessary to release it from the 25(OH)-vitamin D-DBP-complex.

Standards, controls and patient samples which are assayed for 25(OH)-vitamin D are incubated with the releasing reagent. The pre-incubated solutions are then transferred to the microplate coated with 25(OH)-vitamin D, and an anti-25(OH)-vitamin D antibody is added. During an over night incubation step, 25(OH)-vitamin D in the sample and a fixed amount of 25(OH)-vitamin D bound to the microtiter well compete for the binding of the antibody. Then a peroxidase-conjugated antibody is added into each microplate well. A complex of 25(OH)-vitamin D - anti-25(OH)-vitamin D antibody – peroxidase conjugate is formed. Tetramethylbenzidine (TMB) is used as a peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction, whereby the color changes from blue to yellow. The intensity of the yellow color is inversely proportional to the concentration of 25(OH)-vitamin D. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standard. 25(OH)-vitamin D in the samples is determined from this curve.

Test procedure

1. Prior to use in the assay allow all reagents and samples to come to room temperature (18 - 26 °C). For this purpose, open the kit, take out the needed individual components and mix gentle avoiding foam formation
2. Mark the positions of STD (Standards)/SAMPLE/CTRL (Control) on a protocol sheet
3. Label V-tubes (e.g. 1.5 ml Eppendorf-tubes)
4. Pipette **30 µl of STD** (Standard)/**SAMPLE/CTRL** (Control) respectively, into the corresponding tube
5. **Reconstitute RELREAG** (Releasing reagent) (see P. 5., Page 21)
6. Add **300 µl of RELREAG** (releasing reagent) into each tube, vortex shortly

7. **Incubate for 1 hour** at 37 °C in a water bath or heating block (do not use an incubator). Before starting the incubation, ensure that the temperature of the water bath or heating block has reached 37 °C and is maintained constant during the incubation
8. Open tubes carefully and add 600 µL **SAMDIL** (sample dilution buffer). Close the tubes and vortex carefully
9. Take microtiter strips out of the microtiter module. Unused strips must be **covered** with the enclosed foil, stored at 2-8° C and used within 4 weeks
10. Transfer **50 µl of STD** (Standard)/**SAMPLE/CTRL** (Control) from the V-tubes to respective well
11. Add **150 µl of AB** (anti 25(OH)-Vitamin D antibody) into each well
12. Cover the plate tightly with the enclosed foil and **incubate over night** (min. 18 – max. 22 hours) **at 8-10 °C** in the dark
13. Aspirate and wash the wells **5x with 250 µl** of diluted wash buffer. The use of **8-channel pipette** is recommended. Remove remaining wash buffer by hitting the plate against paper towel after the last wash
For TECAN and Dynex instruments a programming protocol can be requested from Immundiagnostik AG
14. Add **200 µl CONJ** (Conjugate) into each well
15. Cover the plate tightly with the enclosed foil and **incubate for 1 hour** at room temperature while shaking
16. Aspirate and wash the wells **5x with 250 µl** of diluted wash buffer. The use of **8-channel pipette** is recommended. Remove remaining wash buffer by hitting the plate against paper towel after the last wash
17. Add **200 µl of SUB** (Substrate) into each well
18. **Incubate for 10 - 15 minutes** at room temperature (18-26°C) in the dark
19. Add **50 µl of STOP** (Stop solution) into each well
20. Determine **absorption** immediately with an ELISA reader at **450 nm**. If the highest extinction of the standards (**STD**) is above the range of the photometer, absorption must be measured immediately at **405 nm** and the obtained results used for evaluation. If possible, the extinctions from each measurement should be compared with extinctions obtained at a reference wavelength, e. g. 595 nm, 620 nm, 630 nm, 650 nm and 690 nm can be used.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend the use of the „4-Parameter-algorithm“.

1. 4-parameter-algorithm

It is recommended a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

9. LIMITATIONS

Samples with 25(OH)-vitamin D concentrations higher than the highest standard should be diluted maximally 1+1 with ready-prepared 1x wash buffer (e. g. 50 µl sample + 50 µl 1x wash buffer) and re-assayed.

Whole blood is not suitable as a sample.

10. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Normal ranges for 25(OH)-Vitamin D₃

Information from ASBMR 2006

Deficiency (seriously deficient)	< 12 ng/ml	resp. < 30 nmol/l
Insufficiency (deficient)	12 - 30 ng/ml	resp. 30 - 75 nmol/l
Sufficiency (adequately supplied)	> 30 ng/ml	resp. > 75 nmol/l

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Reference intervals for 25(OH)-vitamin D₃ (ng/ml)***Males and females**

Age	n	2.5%	97.5%
0 to < 3 months	131	5	42
3 to < 6 months	135	9	60
6 months to < 1 year	147	18	58
1 to < 3 years	394	15	54
3 to < 10 years	619	14	46
10 to < 13 years	286	11	50
13 to < 15 years	275	10	44
15 to < 18 years	390	8	45
> 18 years	421	8	56

*All seasons

(Soldin et al., 2009)

Note

The vitamin D production in the skin is high variable and depends on the season- and daily time, degree of latitude, age, sun protection etc. The normal ranges depend on the method used (e. g. vitamin-D-release from the vitamin D binding protein, DBP) and serve only as orientation.

Literature references

Visser M, Deeg DJ, Puts MT, Seidell JC, Lips P. (2006) Low serum concentrations of 25-hydroxyvitamin D in older persons and the risk of nursing home admission. *Am J Clin Nutr.* Sep;84(3):616-22; quiz 671-2

Grant WB, Holick MF. Benefits and requirements of vitamin D for optimal health: a review. (2005) *Altern Med Rev.* Jun;10(2):94-111. Review

Wicherts IS, van Schoor NM, Boeke AJ, Visser M, Deeg DJ, Smit J, Knol DL, Lips P. (2007) Vitamin D status predicts physical performance and its decline in older persons. *J Clin Endocrinol Metab.* Jun;92(6):2058-65

Soldin OP, Sharma H, Husted L, Soldin SJ. (2009) Pediatric reference intervals for aldosterone, 17alpha-hydroxyprogesterone, dehydroepiandrosterone, testosterone and 25-hydroxy vitamin D3 using tandem mass spectrometry. *Clin Biochem.* Jun;42(9):823-7.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 20)		
Sample	25(OH)-Vitamin D [nmol/L]	CV [%]
1	72.3	7.0

Inter-Assay (n = 20)		
Sample	25(OH)-Vitamin D [nmol/L]	CV [%]
1	84.0	7.0

Specificity

The specificity of the antibody was tested by measuring the cross-reactivity against a range of compounds with structural similarity to 25(OH)-vitamin D₃. The specificity is calculated in per cent, based on the cross-reactivity of these compounds with the anti-25(OH)-vitamin D₃ antibody compared to the 25(OH)-vitamin D₃ antigen:

25(OH)-Vitamin D ₃	100.0 %
25(OH)-Vitamin D ₂	67.8 %
24, 25(OH)-Vitamin D ₃	≥ 100.0 %
Vitamin D ₂ (Ergocalciferol)	0.3 %

Limit of Blank (LoB) = 1.04 ng/ml resp. 2.6 nmol/l

Limit of Detection (LoD) = 1.28 ng/ml resp. 3.2 nmol/l

Limit of Quantitation (LoQ) = 4.80 ng/ml resp. 12 nmol/l

Working range = 4.80 - 96 ng/ml resp. 12 - 240 nmol/l

The evaluation is performed according to the CLSI-Guideline:EP-17-A.

Day-to-Day Variation

Four samples were measured every day within 21 days. The results in nmol/l are shown below:

Sample	Vitamin D mean value [nmol/l]	VC [%]
1	29.8	14.6
2	45.5	13.1
3	70.8	7.3
4	104.3	7.5

Dilution Recovery

To two serum samples were diluted with ready-prepared 1x wash buffer. The samples were analyzed and the results shown below:

Sample	Dilution	Observed [nmol/l]	Expected [nmol/l]	Recovery [%]
A	undiluted	54.5	54.5	
	90%	50.1	49.1	102
	80%	43.4	43.6	100
	70%	36.4	38.2	95
	60%	31.2	32.7	95
	50%	27.4	27.3	101
B	undiluted	77.0	77.0	
	90%	68.0	69.3	98
	80%	59.5	61.6	97
	70%	54.8	53.9	101
	60%	44.6	46.2	97
	50%	31.5	38.5	82



Spiking Recovery

To three serum samples different standard amounts were added. The samples were analyzed and the results shown below:

Sample	Spike [nmol/L]	Target Value [nmol/L]	Obtained Value [nmol/L]	Recovery [%]
6.4	0			
	18.75	25.15	25.80	102.6
	37.50	43.90	35.70	81.3
	75.00	81.40	66.50	81.7
	150.00	156.40	126.90	81.1
	225.00	231.40	204.70	88.5
33.9	0			
	18.75	52.65	49.30	93.6
	37.50	71.40	73.50	102.9
	75.00	108.90	104.80	96.2
	150.00	183.90	185.70	101.0
	225.00	258.90	256.80	99.2
17.2	0			
	18.75	35.95	34.30	95.4
	37.50	54.70	48.00	87.8
	75.00	92.20	81.20	88.1
	150.00	167.20	151.80	90.8
	225.00	242.20	240.80	99.4

Correlation data

Excellent correlation with *direct* ELISA and LC-MS/MS results.

12. PRECAUTIONS

- For *in vitro* diagnostic use only.
- The quality control guidelines should be observed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.

- Reagents of the kit package contain sodium azide and thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic color reactions are described to be also toxic and carcinogenic. Contact with skin or mucous membranes has to be avoided.
- Stop solution consists of sulfuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for in vitro diagnostic use only.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days of receipt of the product. The product shall be send to Immundiagnostik AG together with the complaint in writing.

15. REFERENCES

1. Scharla SH et al. (1996) *Exp Clin Endocrinol Diabetes* 104:289-292
2. Scharla SH et al. (1998) *Osteoporos Int* 8 (Supplement 2):S7-S12
3. Chapuy MC et al. (1996) *J Clin Endocrinol Metab* 81:1129-33
4. Oster P et al. (1983) *Akt Gerontol* 13:221-2
5. Scharla S (1997) Schattauer Verlag, Stuttgart, Seiten 217-242
6. Offermann G (1978) *Dtsch med Wschr* 103:1387-1388
7. Wielders JP, Wijnberg FA. (2009) Preanalytical stability of 25(OH)-vitamin D3 in human blood or serum at room temperature: solid as a rock. *Clin Chem. Aug;55(8):1584-5.*

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

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TECOmedical AG
TECOmedical GmbH