EDI™ Human Fetuin-A ELISA Kit

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LITERATURHINWEISE

- Schinke T, Amendt C, Trindl A, Poschke O, Muller-Esterl W, Jahnen-Dechent W. The serum protein alpha2-HS glycoprotein/Fetuin inhibits apatite formation in vitro and in mineralizing calvaria cells. A possible role in mineralization and calcium homeostasis. *J Biol Chem.* 1996 Aug 23;271:20789-96.
- Jahnen-Dechent W, Schinke T, Trindl A, Muller-Esterl W, Sablitzky F, Kaiser S, Blessing M. Cloning and targeted deletion of the mouse Fetuin gene. *J Biol Chem.* 1997 Dec 12;272:31496-503.
- Mizuno M, Farach-Carson MC, Pinero GJ, Fujisawa R, Brunn JC, Seyer JM, Bousfield GR, Mark MP, Butler WT. Identification of the rat bone 60K acidic glycoprotein as alpha 2HS-glycoprotein. *Bone Miner.* 1991 Apr;13:1-21.
- Price PA, Thomas GR, Pardini AW, Figueira WF, Caputo JM, Williamson MK. Discovery of a high molecular weight complex of calcium, phosphate, Fetuin, and matrix gamma-carboxyglutamic acid protein in the serum of etidronate-treated rats. *J Biol Chem.* 2002 Feb 8;277:3926-34.
- Ketteler M, Vermeer C, Wanner C, Westenfeld R, Jahnen-Dechent W, Floege J. Novel insights into uremic vascular calcification: role of matrix Gla protein and alpha-2-Heremans Schmid glycoprotein/Fetuin. *Blood Purif.* 2002;20:473-476.

- Ketteler M, Bongartz P, Westenfeld R, Wildberger JE, Mahnen AH, Bohm R, Metzger T, Wanner C, Jahnen-Dechent W, Floege J. Association of low Fetuin-A (AHSG) concentrations in serum with cardiovascular mortality in patients on dialysis: a cross-sectional study. *Lancet.* 2003 Mar 8;361:827-833.
- Ketteler M, Wanner C, Metzger T, Bongartz P, Westenfeld R, Gladziwa U, Schurgers LJ, Vermeer C, Jahnen-Dechent W, Floege J. Deficiencies of calcium-regulatory proteins in dialysis patients: a novel concept of cardiovascular calcification in uremia. *Kidney Int Suppl.* 2003 May;S84-87.
- Schäfer C, Heiss A, Schwarz A, Ralf Westenfeld R, Ketteler M, Floege J, Müller-Esterl W, Schinke T, Jahnen-Dechent W. The serum protein alpha 2-Heremans-Schmid glycoprotein/Fetuin-A is a systemically acting inhibitor of ectopic calcification. *J Clin Invest.* 2003 Aug;112(3):357-66.

REFERENZEN

- Stenvinkel P, Wang K, Qureshi AR, Axelsson J, Pecoits-Filho R, Gao P, Barany P, Lindholm B, Jogestrand T, Heimbürger O, Holmes C, Schalling M, Nordfors L. Low fetuin-A levels are associated with cardiovascular death: Impact of variations in the gene encoding fetuin. *Kidney Int.* 2005 Jun;67(6):2383-92.
- Wang AY, Woo J, Lam CW, Wang M, Chan IH, Gao P, Lui SF, Li PK, Sanderson JE. Associations of serum fetuin-A with malnutrition, inflammation, atherosclerosis and valvular calcification syndrome and outcome in peritoneal dialysis patients. *Nephrol Dial Transplant.* 2005 Aug;20(8):1676-85.
- Sato H, Kazama JJ, Wada Y, Kuroda T, Narita I, Gejyo F, Gao P, Yamashita H. Decreased levels of circulating alpha2-Heremans-Schmid glycoprotein/Fetuin-A (AHSG) in patients with rheumatoid arthritis. *Intern Med.* 2007;46(20):1685-91. Epub 2007 Oct 15.
- Lehtinen AB, Burdon KP, Lewis JP, Langefeld CD, Ziegler JT, Rich SS, Register TC, Carr JJ, Freedman BI, Bowden DW. Association of alpha2-Heremans-Schmid glycoprotein polymorphisms with subclinical atherosclerosis. *J Clin Endocrinol Metab.* 2007 Jan;92(1):345-52. Epub 2006 Oct 24
- Lim P, Collet JP, Moutereau S, Guigui N, Mitchell-Heggs L, Loric S, Bernard M, Benhamed S, Montalescot G, Randé JL, Guéret P. Fetuin-A is an independent predictor of death after ST-elevation myocardial infarction. *Clin Chem.* 2007 Oct;53(10):1835-40. Epub 2007 Aug 16
- Ix JH, Wassel CL, Kanaya AM, Vittinghoff E, Johnson KC, Koster A, Cauley JA, Harris TB, Cummings SR, Shlipak MG; Health ABC Study. Fetuin-A and incident diabetes mellitus in older persons. *JAMA.* 2008 Jul 9;300(2):182-8.
- Ix JH, Wassel CL, Bauer DC, Torioan D, Tylavsky FA, Cauley JA, Harris TB, Price PA, Cummings SR, Shlipak MG; for the Health ABC Study. Fetuin-A and Bone Mineral Density in Older Persons: The Health Aging and Body Composition (Health ABC) Study. *J Bone Miner Res.* 2008 Nov 18.
- Giulia Bivona, Chiara Bellia, Antonietta Caruso, Daniela Butera, Bruna Lo Sasso, Patrizia Altavilla, Rosa C. Carollo, Gaia Chiarello and Marcello Ciaccio. Low Serum Fetuin A Levels and Cardiovascular Events in End-Stage Renal Disease (ESRD) Patients. *Research Journal of Medical Sciences* 2 (4): 200-202, 2008
- Mario Cozzolino^a, Andrea Galassi^a, Maria Luisa Biondi^b, Olivia Turri^b, Sergio Papagni^c, Nicola Mongelli^c, Luigi Civita^c, Maurizio Gallieni^a, Diego Brancaccio^a Serum Fetuin-A Levels Link Inflammation and Cardiovascular Calcification in Hemodialysis Patients. *Am J Nephrol* 2006;26:423-429
- G. Metry, P. Stenvinkel, A. R. Qureshi, J. J. Carrero, M. I. Yilmaz, P. Bárány, S. Snaedal, O. Heimbürger, B. Lindholm and M. E. Sulim^{an} Low serum fetuin-A concentration predicts poor outcome only in the presence of inflammation in prevalent haemodialysis patients. *European Journal of Clinical Investigation* 2008;38(11):804 - 811
- Doris Hendig, Veronika Schulz, Marius Arndt, Christiane Szliska, Knut Kleesiek, and Christian Go^o tting1. Role of Serum


- Fetuin-A, a Major Inhibitor of Systemic Calcification, in Pseudoxanthoma Elasticum. Clin Chem 2006;52:227-234
12. G. Marhaug, V. Shah, R. Shroff, H. Varsani, L. R. Wedderburn, C. A. Pilkington and P. A. Brogan. Age-dependent inhibition of ectopic calcification: a possible role for fetuin-A and osteopontin in patients with juvenile dermatomyositis with calcinosis. Rheumatology 2008 47(7):1031-1037
13. Marcello Ciaccio, et al. Changes in serum fetuin-A and inflammatory markers levels in end-stage renal disease (ESRD): effect of a single session haemodialysis. Clin Chem Lab Med 2008;46:212-214
14. J. J. Carrero, P. Stenvinkel, B. Fellström, A. R. Qureshi, K. Lamb, O. Heimbürger, P. Bárány, K. Radhakrishnan, B. Lindholm, I. Soveri, L. Nordfors & P. G. Shiels. Telomere attrition is associated with inflammation, low fetuin-A levels and high mortality in prevalent haemodialysis patients. Journal of Internal Medicine 2007; 263(3):302 – 312
15. C. Fiore, G. Celotta, G. Politi, L. Di Pino, Z. Castelli, R. Mangiafico, S. Signorelli, P. Pennisi. Association of high alpha2-Heremans-Schmid glycoprotein/fetuin concentration in serum and intima-media thickness in patients with atherosclerotic vascular disease and low bone mass. Atherosclerosis 2007;195:110-115
16. Sandro Mazzaferro, Marzia Pasquali, Francesco Pugliese, Giusi Barresi, Iacopo Carbone, Marco Francone, Daniela Sardella, Franco Taggi. Serum Levels of Calcification Inhibition Proteins and Coronary Artery Calcium Score: Comparison between Transplantation and Dialysis. Am J Nephrol 2007;27:75-83.

TECHNICAL ASSISTANCE AND CUSTOMER SERVICE

For technical assistance or place an order, please contact Epitope Diagnostics, Inc. at (858) 693-7877 or fax to (858) 693-7678. www.epitopediagnostics.com

This product is developed and manufactured by
Epitope Diagnostics, Inc.
 San Diego, CA 92126, USA

MDSS
 Burckhardtstrasse 1
 30163 Hannover, Germany

Manufacturer	No. of tests
Catalog Number	Keep away from heat and direct sun light
Concentrate	Store at
In Vitro Diagnostic Device	Use by
Read instructions before use	Lot No.
Authorized Representative In Europe	