

Technical Information

α -GST, Human Urine, Serum, Plasma

Kidney Injury Biomarker for proximal tubular damage

Cat. No.:	TE1056
Tests:	96
Method:	ELISA
Range:	2.5 - 80 $\mu\text{g/l}$
Sensitivity:	1.9 $\mu\text{g/l}$
Incubation time:	2 h
Sample volume:	100 μl
Sample type:	Urine: (Suggested initial dilution 1:2)
Sample preparation:	Overnight or 24 h urine is recommended. As soon as possible after collection sample should be diluted with urine stabilizing buffer (USB). Sample:USB 4:1. After the addition of USB, samples can be stored at 20-25 °C for up to 48 hours, at 2-8 °C for up to one week or at - 20 °C for >1 year. Repeated freeze thawing of samples should be avoided.
	Serum, Plasma: (Suggested initial dilution 1:5) Centrifuge within 3 hours from time of collection and transfer the sample from the original tube for storage. Samples can be stored at 20 – 25 °C for up to 48 hours, at 2 – 8 °C for up to one week or at -20 °C for >1 year. Repeated freeze thawing of samples should be avoided.
Species:	Human
Reference values:	Urine: 0 - 29.0 $\mu\text{g/l}$ Serum, Plasma: 0-12 $\mu\text{g/l}$
Cross reaction:	No cross-reactivity was observed with μ -GST at 500 $\mu\text{g/l}$ or π -GST at 500 $\mu\text{g/l}$

Background:

Alpha Glutathione S-Transferase (α -GST) is a member of the GST superfamily of small cytosolic proteins primarily involved in cellular detoxification reactions. α -GST is approximately 51kDa and can comprise up to 2 % of the cytosolic soluble protein content. As a result of its constitutive nature, α -GST is immediately released into the urine upon lysis of epithelial cells of the proximal tubule.

Intended use:

The α -GST EIA provides a method for the quantitative determination of alpha glutathione S-transferase (α -GST) in human urine. The α -GST EIA can also be used for the determination of α -GST in human, serum and plasma, to assess liver injury.

Elevated urinary α -GST levels are an early indicator of acute kidney injury (AKI) onset in cases of nephrotoxicity, environmental toxicity, cardiothoracic surgery and in transplantation rejection. The release of π -GST has been shown to be associated with distal tubular damage, thus simultaneous measurement of both α -GST and π -GST allows discrimination between proximal and distal tubular damage. Low basal levels of α -GST are released in normal healthy individuals.

For further information please contact / Für weitere Informationen wenden Sie sich bitte an /
Pour plus d'informations, veuillez contacter:

www.tecomedical.com



A EUROBIO SCIENTIFIC COMPANY

Switzerland / Headquarters

TECOmedical AG
Gewerbstrasse 10
4450 Sissach
Phone +41 61 985 81 00
Fax +41 61 985 81 09
Mail info@tecomedical.com

Germany

TECOmedical GmbH
Wasserbreite 57
32257 Bünde
Phone +49 52 23 985 99 99
Fax +49 52 23 985 99 98
Mail info@tecomedical.com

Benelux

TECOmedical Benelux BV
Prins Willem-Alexanderlaan 301
7311 SW Apeldoorn, The Netherlands
Phone +31 30 307 87 30
Fax +31 30 307 49 39
Mail benelux@tecomedical.com

Austria

TECOmedical AG
Phone 0800 20 40 66
Fax 0800 20 40 55
Mail info@tecomedical.com