

## Information

*always your partner*

## EBS SARS-CoV-2 Ag Rapid Test



Lateral flow assay for the quick determination of SARS-CoV-2 antigens in nasopharyngeal swabs

- Ready-to-use extraction units
- 10 minutes to result
- 326 clinically validated samples:
  - Diagnostic sensitivity: 95,7%
  - Diagnostic specificity: 99,1%
- Antigen recognition of the british virus variant




Assay principle	Detection of SARS-CoV-2 Nucleocapsid Antigen
Technology	Qualitative Lateral Flow Immunoassay
Sample material	Nasopharyngeal swabs
Time to result	10 minutes
Analytical sensitivity	Lower Detection Limit: 1,6x10 <sup>2</sup> TCID <sub>50</sub> /ml*
Diagnostic sensitivity	95,7% (110/115) at PCR Cutoff at CT 38
Diagnostic specificity	99,1% (209/211) at PCR Cutoff at CT 38
Description	EBS SARS-CoV-2 Ag Rapid Test
Produced by	EurobioScientific, Les Ulis, France
Packaging	20 tests
Reference	EBS 1020
State Listing	List of COVID-19 antigen tests from the German Federal Institute for Pharmaceuticals and Medical Products (BfArM)
Registration	ANSM (Agence nationale de sécurité de médicament et des produits de santé; France) HAS (Haute Autorité de santé; France)
Test-based approval	Swiss Federal Office of Public Health (BAG)

\* 50% tissue culture infectious dose (TCID<sub>50</sub>): unit for the infectious virus titer

## Assay flow

### 1. SPECIMEN EXTRACTION


- 1.1. Remove the sealing film of the tube with extraction buffer and insert the extraction tube into the tube rack;
- 1.2. Insert the swab into the extraction tube containing the extraction buffer;
- 1.3. Rotate the swab at least 6 times while pressing the tip against the bottom and sides of the extraction tube;
- 1.4. Place the swab in the extraction tube for 1 minute;
- 1.5. Remove the swab from the extraction tube;
- 1.6. Press the nozzle cap tightly onto the tube. The extracted solution will be used as test sample.



### 2. TEST PROCEDURE

- 2.1. Carefully refer to the instruction for use prior to performing the test;
- 2.2. Take out the kits 30 mins before test, ensure that tests and samples are on room temperature;
- 2.3. Place test cassette on flat and clean bench; add 2 drops unknown extracted samples into the sample pad;
- 2.4 Read and record the results after 10 minutes (No longer than 20 minutes). Abnormal results may occur after 20 minutes

### 3. INTERPRETATION of RESULTS



Für weitere Informationen wenden Sie sich bitte an:  
 For further information please contact:  
 Pour plus d'information, veuillez contacter:

[www.tecomedical.com](http://www.tecomedical.com)



A EUROBIO SCIENTIFIC COMPANY

Switzerland / Headquarters  
**TECOmedical AG**  
 Gewerbstrasse 10  
 4450 Sissach  
 Phone +41 61 985 81 00  
 Fax +41 61 985 81 09  
 Mail [info@tecomedical.com](mailto:info@tecomedical.com)

Germany  
**TECOmedical GmbH**  
 Wasserbreite 57  
 32257 Bünde  
 Phone +49 52 23 985 99 99  
 Fax +49 52 23 985 99 98  
 Mail [info@tecomedical.com](mailto:info@tecomedical.com)

Benelux  
**TECOmedical Benelux BV**  
 Prins Willem-Alexanderlaan 301  
 7311 SW Apeldoorn, The Netherlands  
 Phone +31 30 307 87 30  
 Fax +31 30 307 49 39  
 Mail [benelux@tecomedical.com](mailto:benelux@tecomedical.com)

Austria  
**TECOmedical AG**  
 Phone 0800 20 40 66  
 Fax 0800 20 40 55  
 Mail [info@tecomedical.com](mailto:info@tecomedical.com)