

**Certificate of Analysis**

**Analysezertifikat**

**Certificat d'Analyse**

<b>Product name</b>	Complement factor B/Ba, Human, clone M20/6		
<b>Catalog number</b>	TE2255-100UG		
<b>Lot number</b>	xxxxxXxxxx	<b>Expiry date</b>	MMM YYYY
<b>Volume</b>	1 ml	<b>Amount</b>	100 µg
<b>Formulation</b>	0.2 µm filtered in PBS+0.1%BSA+0.02%NaN3	<b>Concentration</b>	100 µg/ml
<b>Host Species</b>	Mouse IgG1	<b>Conjugate</b>	None
<b>Endotoxin</b>	N.A.	<b>Purification</b>	Protein G
<b>Storage</b>	4°C		

**Application notes**

	IHC-F	IHC-P	IF	FC	FS	IA	IP	W
Reference #								
Yes						•	•	•
No								
N.D.	•	•	•	•	•			

*N.D.= Not Determined; IHC = Immuno histochemistry; F = Frozen sections; P = Paraffin sections; IF = Immuno Fluorescence; FC = Flow Cytometry; FS = Functional Studies; IA = Immuno Assays; IP = Immuno Precipitation; W = Western blot*

Dilutions to be used depend on detection system applied. It is recommended that users test the reagent and determine their own optimal dilutions. The typical starting working dilution is 1:50.

**General Information**

<b>Description</b>	The monoclonal antibody M20/6 recognizes human complement Factor B/Ba. Factor B is an acute-phase protein. Levels of factor B increase during inflammation. Complement factor B is a single- chain molecule of 764 amino acids (MW of 90 kD), including a leader peptide of 25 amino acids. Factor B provides the catalytic subunit of the C3/C5 convertases of the alternative complement pathway. Assembly of the C3 convertase (C3bBb) requires binding of factor B to C3b (C3.H2O) and factor D-mediated cleavage of bound factor B resulting in the release of Ba (MW 33 kD). The C3 convertase is stabilized by the binding of properdin. This provides a positive amplification loop for the classical and alternative complement pathways. Bb (MW 60 kD) is the serine protease element of this convertase. After cleavage of C3, the C5 convertase ((C3b)2Bb) is formed. The Bb fragment may be regarded as a better indicator of the alternative pathway of complement activation than Ba as impaired renal filtration does directly influence Ba levels. Whereas Bb is elevated 2.2-fold in chronic renal failure (CRF) patients, plasma concentrations of Ba are 8.4-fold higher in CRF and 16-fold higher in end-stage renal disease (ESRD) patients compared to normals. In addition to complement activation, factor B fragments participate in other immunological functions.
<b>References</b>	<ol style="list-style-type: none"> <li>1. Oppermann, M et al; Complement in patients with renal failure as detected through the quantitation of fragments of the complement proteins C3, C5 and factor B. <i>Klin Wochenschr</i> 1988, 66: 857</li> <li>2. Oppermann M. et al; Quantitation of components of the alternative pathway of complement (APC) by enzyme-linked immunosorbent assays. <i>J Immunol Methods</i> 1990, 133: 181</li> <li>3. Oppermann M. et al; Elevated plasma levels of the immunosuppressive complement fragment Ba in renal failure. <i>Kidney Int</i> 1991, 40: 939</li> </ol>
<b>Storage&amp;stability</b>	Product should be stored at 4°C. Under recommended storage conditions, product is stable for at least one year.
<b>Precautions</b>	For research use only. Not for use in or on humans or animals or for diagnostics. It is the responsibility of the user to comply with all local/state and federal rules in the use of this product. Hycult Biotech is not responsible for any patent infringements that might result from the use or derivation of this product.

---

We hereby certify that the above-stated information is correct and that this product has been successfully tested by the Quality Control Department. This product was released for sale according to the existing specifications. This document has been produced electronically and is valid without a signature.

Approved by Manager of QC

Date

Do you have any questions or comments regarding this product? Please contact us via [info@tecomedical.com](mailto:info@tecomedical.com)

[www.tecomedical.com](http://www.tecomedical.com)



A EURO BIO SCIENTIFIC COMPANY

Switzerland  
**TECO medical AG**  
Gewerbstrasse 10  
4450 Sissach  
Phone +41 61 985 81 00  
Fax +41 61 985 81 09  
Mail [info@tecomedical.com](mailto:info@tecomedical.com)

Germany  
**TECO medical GmbH**  
Wasserbreite 57  
32257 Bünde  
Phone +49 52 23 985 99 99  
Fax +49 52 23 985 99 98  
Mail [info@tecomedical.com](mailto:info@tecomedical.com)

Benelux  
**TECO medical Benelux BV**  
Computerweg 22  
3542 DR Utrecht, The Netherlands  
Phone +31 30 307 87 30  
Fax +31 30 307 49 39  
Mail [benelux@tecomedical.com](mailto:benelux@tecomedical.com)

Austria  
**TECO medical AG**  
Phone 0800 20 40 66  
Fax 0800 20 40 55  
Mail [info@tecomedical.com](mailto:info@tecomedical.com)