

**Certificate of Analysis**

**Analysezertifikat**

**Certificat d'Analyse**

<b>Product name</b>	C6, Human, clone WU 6-4		
<b>Catalog number</b>	TE2276-100UG		
<b>Lot number</b>	xxxxxXxxxx	<b>Expiry date</b>	MMM YYYY
<b>Volume</b>	1 ml	<b>Amount</b>	100 µg
<b>Formulation</b>	0.2 µm filtered in PBS+0.1%BSA	<b>Concentration</b>	100 µg/ml
<b>Host Species</b>	Mouse IgG1	<b>Conjugate</b>	None
<b>Endotoxin</b>	<24 EU/mg	<b>Purification</b>	Protein G
<b>Storage</b>	4°C		

**Application notes**

	IHC-F	IHC-P	IF	FC	FS	IA	IP	W
Reference #								
Yes		•	•	•	•	•	•	•
No								
N.D.	•							

*N.D.= Not Determined; IHC = Immuno histochemistry; F = Frozen sections; P = Paraffin sections; IF = Immuno Fluorescence; FC = Flow Cytometry; FS = Functional Studies; IA = Immuno Assays; IP = Immuno Precipitation; W = Western blot*

Dilutions to be used depend on detection system applied. It is recommended that users test the reagent and determine their own optimal dilutions. The typical starting working dilution is 1:50. For functional studies, in vitro dilutions have to be optimized in user's experimental setting.

**General Information**

<b>Description</b>	The monoclonal antibody WU 6-4 recognizes C6, one of the components of the terminal complement complex (TCC), also known as the membrane attack complex (MAC). Proteolytic cleavage of C5 by C5 convertase generates C5b which initiates assembly of the C5b-9 MAC. This complex is assembled from five precursor molecules in the serum, finalized with the polymerization of C9 which accompanies insertion of the complex into the cell membrane causing cellular lysis. Polymorphonuclear leukocytes represent a major source of C6. The monoclonal antibody WU 6-4 recognizes native C6 but not C6 after incorporation into the TCC. Furthermore, WU 6-4 is capable of inhibiting C5b6 and TCC formation in the fluid phase and to a lesser extent, haemolysis. Haemolytic activity could be inhibited only when WU 6-4 was introduced before C5b6 formation.
<b>Immunogen</b>	Human C6.
<b>Cross reactivity</b>	Primate: Yes; Rat: Yes.
<b>References</b>	<ol style="list-style-type: none"> <li>Würzner R et al; Inhibition of terminal complement complex formation and cell lysis by monoclonal antibodies. <i>Complement and Inflammation</i> 1991, 8: 328</li> <li>Würzner R et al; Functionally active complement proteins C6 and C7 detected in C6- or C7- deficient individuals. <i>Clin and Exp Immunol</i> 1991, 83: 430</li> <li>Würzner R et al; Blood dendritic cells carry terminal complement complexes on their cell surface as detected by new developed neoepitope-specific monoclonal antibodies. <i>Immunology</i> 1991, 174:32</li> <li>Orren A et al; Properties of a low molecular weight complement component C6 found in human subjects with subtotal C6 deficiency. <i>Immunol</i> 1992, 75: 10</li> <li>Würzner R et al; Importance of the third thrombospondin repeat of C6 for terminal complement complex assembly. <i>Immunol</i> 1995, 85: 214</li> <li>Würzner R et al. (1995); Molecular basis of subtotal complement C6 deficiency: a carboxyterminally truncated but functionally active C6. <i>J of Clin Invest</i> 1995, 95: 1877</li> </ol>
<b>Storage&amp;stability</b>	Product should be stored at 4°C. Under recommended storage conditions, product is stable for at least one year.

**Precautions** For research use only. Not for use in or on humans or animals or for diagnostics. It is the responsibility of the user to comply with all local/state and federal rules in the use of this product. Hycult Biotech is not responsible for any patent infringements that might result from the use or derivation of this product.

---

We hereby certify that the above-stated information is correct and that this product has been successfully tested by the Quality Control Department. This product was released for sale according to the existing specifications. This document has been produced electronically and is valid without a signature.

Approved by Manager of QC

Date

Do you have any questions or comments regarding this product? Please contact us via [info@tecomedical.com](mailto:info@tecomedical.com)

[www.tecomedical.com](http://www.tecomedical.com)



A EURO BIO SCIENTIFIC COMPANY

Switzerland  
**TECO medical AG**  
Gewerbstrasse 10  
4450 Sissach  
Phone +41 61 985 81 00  
Fax +41 61 985 81 09  
Mail [info@tecomedical.com](mailto:info@tecomedical.com)

Germany  
**TECO medical GmbH**  
Wasserbreite 57  
32257 Bünde  
Phone +49 52 23 985 99 99  
Fax +49 52 23 985 99 98  
Mail [info@tecomedical.com](mailto:info@tecomedical.com)

Benelux  
**TECO medical Benelux BV**  
Computerweg 22  
3542 DR Utrecht, The Netherlands  
Phone +31 30 307 87 30  
Fax +31 30 307 49 39  
Mail [benelux@tecomedical.com](mailto:benelux@tecomedical.com)

Austria  
**TECO medical AG**  
Phone 0800 20 40 66  
Fax 0800 20 40 55  
Mail [info@tecomedical.com](mailto:info@tecomedical.com)