

Certificate of Analysis

Analysezertifikat

Certificat d'Analyse

Product name	C9 (C-terminal), Human, clone X197		
Catalog number	TE2111-100UG		
Lot number	xxxxxXxxxx	Expiry date	MMM YYYY
Volume	1 ml	Amount	100 µg
Formulation	0.2 µm filtered in PBS+0.1%BSA	Concentration	100 µg/ml
Host Species	Mouse IgG1	Conjugate	None
Endotoxin	<24 EU/mg	Purification	Protein G
Storage	4°C		

Application notes

	IHC-F	IHC-P	IF	FC	FS	IA	IP	W
Reference #								
Yes	•	•		•	•	•		•
No								
N.D.			•				•	

N.D.= Not Determined; IHC = Immuno histochemistry; F = Frozen sections; P = Paraffin sections; IF = Immuno Fluorescence; FC = Flow Cytometry; FS = Functional Studies; IA = Immuno Assays; IP = Immuno Precipitation; W = Western blot

Dilutions to be used depend on detection system applied. It is recommended that users test the reagent and determine their own optimal dilutions. The typical starting working dilution is 1:50.

- FS: inhibition of C9 binding to C5b-8.

General Information

Description	Human C9 is a soluble glycoprotein of 61 kD. It is the last component in the assembly of the membrane attack complex (MAC). When the complement is activated on target membranes, multiple copies of C9 bind to the C5b-8 complex and assemble the barrel-shaped pore which causes cell lysis. The conformation of C9 changes from globular to a tubular form. The binding of C9 to C5b-8 plays a key role in the function of C9.
References	1. Hatanaka, M et al; Analysis of C5b-8 binding sites in the C9 molecule using monoclonal antibodies: participation of two separate epitopes of C9 in C5b-8 binding. Mol Immunol 1992, 29: 911
Storage&stability	Product should be stored at 4°C. Under recommended storage conditions, product is stable for at least one year.
Precautions	For research use only. Not for use in or on humans or animals or for diagnostics. It is the responsibility of the user to comply with all local/state and federal rules in the use of this product. Hycult Biotech is not responsible for any patent infringements that might result from the use or derivation of this product.

We hereby certify that the above-stated information is correct and that this product has been successfully tested by the Quality Control Department. This product was released for sale according to the existing specifications. This document has been produced electronically and is valid without a signature.

Approved by Manager of QC

Date

Do you have any questions or comments regarding this product? Please contact us via info@tecomedical.com

www.tecomedical.com



A EUROBIO SCIENTIFIC COMPANY

Switzerland
TECOmedical AG
Gewerbestrasse 10
4450 Sissach
Phone +41 61 985 81 00
Fax +41 61 985 81 09
Mail info@tecomedical.com

Germany
TECOmedical GmbH
Wasserbreite 57
32257 Bünde
Phone +49 52 23 985 99 99
Fax +49 52 23 985 99 98
Mail info@tecomedical.com

Benelux
TECOmedical Benelux BV
Computerweg 22
3542 DR Utrecht, The Netherlands
Phone +31 30 307 87 30
Fax +31 30 307 49 39
Mail benelux@tecomedical.com

Austria
TECOmedical AG
Phone 0800 20 40 66
Fax 0800 20 40 55
Mail info@tecomedical.com