

Technical Information

Sclerostin TECO[®] High sensitive Regulation of Bone Turnover

Cat. No.: TE1023-HS
 Tests: 96
 Method: **ELISA**
 Range: 0.05 – 3 ng/ml
 Sensitivity: LOD: 0.008 ng/ml
 LLOQ: 0.01 ng/ml
 ULOQ: 3.34 ng/ml
 Incubation time: 4.5 hours
 Sample volume: 25 µl
 Sample type: Plasma, Serum, Cell culture
 Sample preparation: Samples should be collected without hemolysis. Allow the blood to clot, separate serum by centrifugation and test samples within 2 hours. Serum is stable for 3 days at room temperature, 5 days at 2 – 8 °C, 2 years at -20 °C. Longer storage at -80°C or colder. Maximum 3 freeze- and thaw cycles.

Species: Human

Reference values:

Subjects	Mean (ng/ml)	SD (ng/ml)	Mean (pmol/l)	SD (pmol/l)	N
Pre-menopausal female	0.48	0.14	21.12	6.16	20
Post-menopausal female	0.58	0.17	25.52	7.48	19
Men	0.64	0.15	38.16	6.6	10

Clinically Healthy Subjects

Intended use

- Sclerosteosis
- Van Buchem disease
- Osteoporosis
- Chronic Kidney Disease
- Cardiac Research
- Others

For further information please contact / Für weitere Informationen wenden Sie sich bitte an /
Pour plus d'informations, veuillez contacter:

www.tecomedical.com



A EUROBIO SCIENTIFIC COMPANY

Switzerland / Headquarters

TECO medical AG
Gewerbstrasse 10
4450 Sissach
Phone +41 61 985 81 00
Fax +41 61 985 81 09
Mail info@tecomedical.com

Germany

TECO medical GmbH
Wasserbreite 57
32257 Bünde
Phone +49 52 23 985 99 99
Fax +49 52 23 985 99 98
Mail info@tecomedical.com

Benelux

TECO medical Benelux BV
Prins Willem-Alexanderlaan 301
7311 SW Apeldoorn, The Netherlands
Phone +31 30 307 87 30
Fax +31 30 307 49 39
Mail benelux@tecomedical.com

Austria

TECO medical AG
Phone 0800 20 40 66
Fax 0800 20 40 55
Mail info@tecomedical.com