

Information

always your partner

TECO[®] Fast Candida Ag & IgG



2 lateral flow assays for rapid, parallel diagnosis of invasive candidiasis and candidemia on the antigen and antibody level

- Fluorescence-based lateral flow assays designed for readout via reader
- (Semi-)quantitative, calculation of IgG concentration in AU/mL and Mannan Antigen in pg/mL
- Validated for serum samples
- 2 results in 40 minutes obtained from 305 μ L of serum
- Detection limit IgG: 7 AU/mL, cut-off value of 135 AU/mL
- Detection limit Antigen: 21.5 pg/mL, cut-off value of 60 pg/mL
- Increase in diagnostic significance through combination of two parameters

	Fast Candida Mannan IgG	Fast Candida Mannan Antigen
Clinical sensitivity	86.2% (n = 222)	88.1% (n = 155)
Clinical specificity	95.3% (n = 222)	96.9% (n = 155)



TECO® *Fast* Candida IgG and Antigen LFA – sample preparation



TECO® *Fast* Candida Mannan Antigen LFA
Catalogue No.: TE1081

&

TECO® *Fast* Candida IgG Antibody LFA
Catalogue No.: TE1083

Für weitere Informationen wenden Sie sich bitte an:
For further information please contact:
Pour plus d'information, veuillez contacter:

www.tecomedical.com



A EUROBIO SCIENTIFIC COMPANY

Switzerland / Headquarters

TECO medical AG
Gewerbstrasse 10
4450 Sissach
Phone +41 61 985 81 00
Fax +41 61 985 81 09
Mail info@tecomedical.com

Germany

TECO medical GmbH
Wasserbreite 57
32257 Bünde
Phone +49 52 23 985 99 99
Fax +49 52 23 985 99 98
Mail info@tecomedical.com

Benelux

TECO medical Benelux BV
Prins Willem-Alexanderlaan 301
7311 SW Apeldoorn, The Netherlands
Phone +31 30 307 87 30
Fax +31 30 307 49 39
Mail benelux@tecomedical.com

Austria

TECO medical AG
Phone 0800 20 40 66
Fax 0800 20 40 55
Mail info@tecomedical.com