

## Information

*always your partner*

### 2019-nCOV IgG/IgM Rapid Test




Lateral flow assay for the qualitative determination of 2019-nCOV IgG/IgM

- Combo-Design: IgM/IgG determination on one test cassette;
- Flexible usage of serum, plasma and whole blood samples;
- Turn-Around Time of 10 minutes;
- Excellent performance confirmed in multiples international studies with more than 1`600 samples

Test principle	Determination of IgM and IgG against SARS-CoV-2 Nucleocapsid and Spike protein
Technology	Lateral flow immunoassay on test cassettes
Sample type	Serum, plasma and whole blood
Incubation time	10 minutes
Positive Predicted Value (PPV)	91.5% (n=165)
Negative Predicted Value (NPV)	96.3% (n=297)
Sensitivity	93.2% (n=162)
Specificity	95.5% (n=300)
Description	2019-nCoV IgG/IgM Rapid Test
Produced by	Dynamiker Biotechnology Co., Ltd.; Tianjin, China
Packaging	40 test cassettes
Reference	DNK-1419-1

## Assay Flow



Für weitere Informationen wenden Sie sich bitte an:  
 For further information please contact:  
 Pour plus d'information, veuillez contacter:

[www.tecomedical.com](http://www.tecomedical.com)



A EUROBIO SCIENTIFIC COMPANY

Switzerland / Headquarters

**TECOmedical AG**

Gewerbestrasse 10

4450 Sissach

Phone +41 61 985 81 00

Fax +41 61 985 81 09

Mail [info@tecomedical.com](mailto:info@tecomedical.com)

Germany

**TECOmedical GmbH**

Wasserbreite 57

32257 Bünde

Phone +49 52 23 985 99 99

Fax +49 52 23 985 99 98

Mail [info@tecomedical.com](mailto:info@tecomedical.com)

Benelux

**TECOmedical Benelux BV**

Computerweg 22

3542 DR Utrecht, The Netherlands

Phone +31 30 307 87 30

Fax +31 30 307 49 39

Mail [benelux@tecomedical.com](mailto:benelux@tecomedical.com)

Austria

**TECOmedical AG**

Phone 0800 20 40 66

Fax 0800 20 40 55

Mail [info@tecomedical.com](mailto:info@tecomedical.com)