

## Information

*always your partner*

# EDI™ Novel Coronavirus COVID-19 IVD CE Nucleocapsid IgG ELISA Kit

ELISA test for the qualitative determination of COVID-19 IgG in serum samples

- Microtiterplates coated with recombinant COVID-19 full length nucleocapsid protein;
- Positive predicted value of 98.4 %;
- Negative predicted value of 99.8%
- Very high specificity of 99.8%

Test principle	Qualitative determination of human IgG's against COVID-19
Technology	Indirect Immunoassay using microtiterplates coated with COVID-19 recombinant full length nucleocapsid protein
Sample type	Serum
Incubation time	80 minutes
Positive Predicted Value (PPV)	98.4% (n=187)
Negative Predicted Value (NPV)	99.8% (n=624)
Specificity	99.8% (n=570)
Description	EDI™ Novel Coronavirus COVID-19 IgG ELISA Kit
Produced by	Epitope Diagnostics, Inc, San Diego, USA
Packaging	1 Microtiterplate with 96 wells
Reference	KT-1032

Für weitere Informationen wenden Sie sich bitte an:  
For further information please contact:  
Pour plus d'information, veuillez contacter:

[www.tecomedical.com](http://www.tecomedical.com)



---

**A EUROBIO SCIENTIFIC COMPANY**

**Switzerland / Headquarters**

**TECOmedical AG**

Gewerbestrasse 10

4450 Sissach

Phone +41 61 985 81 00

Fax +41 61 985 81 09

Mail [info@tecomedical.com](mailto:info@tecomedical.com)

**Germany**

**TECOmedical GmbH**

Wasserbreite 57

32257 Bünde

Phone +49 52 23 985 99 99

Fax +49 52 23 985 99 98

Mail [info@tecomedical.com](mailto:info@tecomedical.com)

**Benelux**

**TECOmedical Benelux BV**

Computerweg 22

3542 DR Utrecht, The Netherlands

Phone +31 30 307 87 30

Fax +31 30 307 49 39

Mail [benelux@tecomedical.com](mailto:benelux@tecomedical.com)

**Austria**

**TECOmedical AG**

Phone 0800 20 40 66

Fax 0800 20 40 55

Mail [info@tecomedical.com](mailto:info@tecomedical.com)