

Information

always your partner

EurobioPlex FluCoSyn



Differential diagnosis of respiratory infections
(Influenza A/B, COVID-19, RSV A/B) by RT-PCR

- Manufactured in France
- Screening of 3 viruses in 1 well
- Very high sensitivity and specificity
- Validated on upper and lower respiratory tract samples

Test principle	Influenza A/B, SARS-CoV-2 and RSV A/B Pathogen Detection
Technology	Real-time RT-PCR
Targeted genes	RdRp gene, N gene plus genomic fragments of RSV and influenza
Sample type	Nasopharyngeal swab, Nasopharyngeal aspirate, Bronchial aspiration, Bronchoalveolar fluid, Sputum
Controls included	Positive control, negative control and internal control
Duration process	1 hour 30 minutes
Validated on	LC480 [®] , CFX96 [®] , ABI7500
Clin. Sensitivity	≥99% (Flu A/B = 97%)
Clin. Specificity	>99%
Positive Predicted Value	>99%
Negative Predicted Value	≥98% (Flu A/B = 96%)
Description	EurobioPlex FluCoSyn Multiplex
Produced by	Eurobio-Scientific, Les Ulis, France
Packaging	192 reactions
Reference	EBX-042-192

Für weitere Informationen wenden Sie sich bitte an:
For further information please contact:
Pour plus d'information, veuillez contacter:

www.tecomedical.com



A EUROBIO SCIENTIFIC COMPANY

Switzerland / Headquarters

TECOmedical AG

Gewerbestrasse 10

4450 Sissach

Phone +41 61 985 81 00

Fax +41 61 985 81 09

Mail info@tecomedical.com

Germany

TECOmedical GmbH

Wasserbreite 57

32257 Bünde

Phone +49 52 23 985 99 99

Fax +49 52 23 985 99 98

Mail info@tecomedical.com

Benelux

TECOmedical Benelux BV

Computerweg 22

3542 DR Utrecht, The Netherlands

Phone +31 30 307 87 30

Fax +31 30 307 49 39

Mail benelux@tecomedical.com

Austria

TECOmedical AG

Phone 0800 20 40 66

Fax 0800 20 40 55

Mail info@tecomedical.com